Automatic Control Systems

**B-ACS 60/120
Instruction Manual**

Version 002

BRUKER

The information in this manual may be altered without notice.

BRUKER accepts no responsibility for actions taken as a result of use of this manual. BRUKER accepts no liability for any mistakes contained in the manual, leading to coincidental damage, whether during installation or operation of the instrument. Unauthorised reproduction of manual contents, without written permission from the publishers, or translation into another language, either in full or in part, is forbidden.

This manual was written by

Stanley J. Niles

© February 4, 1997: Bruker Elektronik GmbH

Rheinstetten, Germany

P/N: Z31285
DWG-Nr: 997002

DECLARATION OF CONFORMITY

The undermentioned product

B-ACS 60 Sample Changer H1080

conforms to the main requirements set by the commission for the Harmonization of Regulations of the EU Member States with regards to electromagnetic compatibility (EMI: 89/336/EWG) and safety (Low Voltage Electrical Equipment: 73/23/EWG) regulations.

For the assessment the following norms were applied:

EMI: EN 55 011; EN 50 082-1

Safety: EN 61 010-1

Test report UNI KA: 9522
Documentation: Z31285 Docu Standard: Sample Changer


Manufacturer's Name: BRUKER ELEKTRONIK
Manufacturer's Address: 76287 Rheinstetten, Akazienweg 2, Germany

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(Signature)

Contents

	Contents	ix
1	Introduction (ENGLISH)	13
1.1	The BRUKER Automatic Sample Changer	13
1.2	Site Considerations	13
1.3	Contact for Additional Technical Assistance	13
2	Safety Considerations	15
2.1	Before Mounting the Automatic Sample Changer	15
2.2	While Mounting the Column	15
2.3	During Operation	15
2.4	Potentially Hazardous Areas	15
2.5	Extreme Temperatures	16
3	Mounting Instructions	17
3.1	Caution	17
3.2	Mounting the Light Barrier Assembly	18
3.3	Mounting the Cabinet	18
3.4	Preparing the Column for Assembly	19
3.5	Assembling the Column	20
3.6	Mounting the Column Assembly to the Magnet	21
3.7	Mounting the Pneumatic Arm Assembly	23
3.8	Securing the Bottom of the Column Assembly	24
3.9	Connecting the Hoses and Cables	25
3.10	Description of Input and Output Locations	27
4	Special Tools	29
4.1	Special Tools for the Automatic Sample Changer	29
5	Control Elements	31
5.1	General	31
5.2	Error Messages	31
5.3	Errors Requiring User Intervention	31
5.4	Communication Errors	33
6	Operator Maintenance	35
6.1	Greasing the Horizontal Arm	35
6.2	Greasing the Vertical Cylinder	35
6.3	Arm Adjustment	35
6.4	Cabinet Adjustment	35
6.5	Vertical Cylinder	35

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Introduction

1

The BRUKER Automatic Sample Changer

1.1

The BRUKER Automatic Sample Changer, used in conjunction with BRUKER DISNMR, UXNMR, or XWINNMR software, provides dialog-guided facilities which allow the user to easily and effectively perform automatic (continuous) experiments. Features include a 60 or 120 sample capacity, random accessing of samples, positive sample identification with the optional barcode reader, and temperature control of individual samples with the optional sample heater unit.

The standard NMR software includes the comprehensive automation package required for the Automatic Sample Changer. It features automatic instrument optimization, phasing, integration and plotting. Specialized routines assist data interpretation. Easy set-up procedures are accomplished via dialog software at three user-interface levels.

Site Considerations

1.2

The Automatic Sample Changer should be setup in a standard laboratory environment. Maximum room temperature should not exceed 30°C. For more information refer to the Avance Spectrometer manual on site planning.

Contact for Additional Technical Assistance

1.3

For further technical assistance, please contact us at:

BRUKER Analytische Messtechnik GmbH
Am Siberstreifen
D-76287 Rheinstetten
Tel.: [+49] 721 5161 0
Fax: [+49] 721 5171 01

Introduction

E
N
G
L
I
S
H

Safety Considerations

2

Before Mounting the Automatic Sample Changer

2.1

Make sure the magnet is firmly secured to its base in order to prevent the instrumentation from tipping over.

While Mounting the Column

2.2

When assembling the column, do not tighten the top fastening ring too hard, as this may cause damage to the column surfaces.

When fixing the column to the magnet, mount both bottom arms to prevent obstruction when manipulating or changing probeheads.

When handling tools, screws, or any metallic parts beware of the strong magnetic field.

During Operation

2.3

Excessive tension of the magazine belt will damage the belt and wheel bearings. (Refer to the operating instructions in the Automatic Sample Changer manual for adjustment of the tension of the magazine belt).

Beware of the strong magnetic field while working around the magnet. Keep all metal objects, such as tools, screws, or any metallic parts away from the magnet. Remove any mechanical watches or metallic objects while working around the magnet.

When the Automatic Sample Changer is running, avoid putting hands or objects in the path of the arm, magazine belt or pinchers, as this may cause personal injury or damage to the equipment.

Potentially Hazardous Areas

2.4

The symbol shown on the left indicates a **potentially hazardous area** (ISO 3864; DIN 40008).

The symbol is placed on the following areas of the Automatic Sample Changer:

On the upper left half of the front side (side that is away from the magnet) of the vertical cylinder. **Warning:** When the sample changer is in operation keep hands and other objects away from the pinchers on the end of the vertical cylinder.


1. On the front side of the upper column assembly. **Warning:** Keep hands and other objects away from the path of the pneumatic arm and magazine belt (figure 3.7) during operation of the sample changer.

Extreme Temperatures

2.5

The sample warmer located on the magazine belt reaches extreme temperatures: **Warning:** Do not insert your fingers into the sample warmer. When the sample warmer is operating it can reach temperatures exceeding 70°C (figure 2.1).

Figure 2.1. Extreme Temperature Warning:

Warning!! Extreme Temperature! Do not put your finger into the sample warmer!

Mounting Instructions


3

Caution

3.1

Before mounting the Automatic Sample Changer, it is important to make sure that the magnet is firmly anchored to its base in order to prevent the instrumentation from tipping over. When doing so, carefully align the legs and lifting lugs used to secure the Sample Changer. The Sample Changer is held to the magnet by four adjustable arm assemblies. The arm assemblies should be aligned with the magnet as shown in figure 3.1.

Figure 3.1. Location of Adjustable Arm Assemblies


Mounting the Light Barrier Assembly

3.2

Note: Magnets with the BOSS 1 and BOSS 2 Shim Systems do not require this step, they are delivered with the light barrier cylinder already installed. Mount the light barrier cylinder by placing it onto the top of the shim system and fastening the three side screws (figure 3.2).

Figure 3.2. Light Barrier for the Shim System

Mounting the Cabinet

3.3


The sample changer cabinet is fitted with a display unit on the front side. The rear side is the side that faces the magnet. Remove the two rear side panels from the rear of the cabinet. Loosen the screws from the top cylindrical collar (figure 3.3) before mounting the unit onto the column.

Place the cabinet on its side on a soft surface such as a blanket or carpet. This position will make it easier to slide the column through the cabinet. Open the front panel of the cabinet by loosening the screws.

Note: The top end of the round column piece is the end with the shortest distance from the side opening to the end of the column piece.

Insert the top end of the round column piece from the bottom side of the cabinet into the cabinet cylindrical collar, pushing it through until the side opening on the column disappears inside the cabinet (figure 3.3). Mount the cabinet resting ring between the side opening and the cylindrical collar. Tighten both the resting ring and the collar in the position shown in figure 3.3.

Figure 3.3. Front view of the Automatic Sample Changer

Preparing the Column for Assembly

3.4

Slide the clamp cover and the square clamp (figure 3.4) over the bottom of the round column towards the bottom of the cabinet. The clamp cover should be next to the cabinet and the square clamp underneath the cover. Do not tighten the clamp.

To prepare the square column piece for mounting, pull the cables and pneumatic hoses out of the inside of the top of the square column piece. Insert these cables and hoses (upwards) through the round column piece. Reach through the cabinet opening and the side opening of the round column piece and pull the cables and hoses through until the ends are outside the side opening of the round piece (figure 3.4).

Measure the distance from the top of the shim system of the magnet to the ground. If the Column is fitted with a base plate, then measure from the top of the shim system to the top of the base plate (figure 3.1). In either case subtract 0,5 cm from the measurement. This is the length that the column needs to be from the bottom of the square column piece to the top of the belt of the cabinet. You will need this measurement for the assembly of the column.

Assembling the Column

3.5


Refer to figure 3.4 for this step.

Assemble the two column pieces, placing the square piece over the round piece. Using the measurement obtained in section 3.4 align the distance of the bottom of the square piece to the top of the belt of the cabinet. It is **important** for the proper operation of the sample changer that this length matches the distance from the ground (or the top of the base plate for the column assembly) to the top of the shim system minus 0,5 cm.

Fasten the two column pieces together by placing the four vertical locking mechanisms in the inside corners of the square column piece. Recheck the length of the column as described above, and adjust if necessary.

Slide the square clamp downwards until it sets over the locking mechanisms. Tighten the clamp securely. Slide the clamp cover downwards until it sets firmly against the square clamp.

Figure 3.4. Column Assembly

Mounting the Column Assembly to the Magnet

3.6


Fasten the two top adjustable arm assemblies to the round portion of the column. Remove the bolts and brackets (figure 3.5) from the magnet end of the adjustable arms and place them next to the magnet lifting lugs.

Lift the column assembly upwards so that the cabinet is on the top. Secure the two top adjustable arms to the lifting lugs of the magnet with the bolts and brackets. The column should be turned in such a way that the connections at the bottom are directed towards the magnet (figure 3.1). **Tip:** It is easier to accomplish this step with two people, one to hold the column vertical and the other to secure the adjustable arms.

Mounting Instructions


Note: If the magnet has Vibration Dampers and the column uses a support pictured in figure 3.1, then the column must be leveled at this time (as described in section 3.8) before mounting the pneumatic arm assembly.

Figure 3.5. Adjustable Arm Assembly

E
N
G
L
I
S
H

Figure 3.6. Arrangement of Arm and Cabinet Assemblies

Mounting the Pneumatic Arm Assembly

3.7


Slide the pneumatic arm resting ring (figure 3.6) over the top of the round column piece and let it rest on the top of the cabinet (don't tighten it).

Place the pneumatic arm assembly over the top of the column, sliding the two fixture rings (figure 3.6) over the round column piece. Tighten the rings just enough to keep the pneumatic arm assembly from sliding down.

Rotate the pneumatic arm assembly (figure 3.6) until the pincher is directly over the shim system (figure 3.7). The distance between the bottom of the pincher (when the cylinder is in the down position) and the top of the shim system should be 3 cm. Adjust this distance by loosening the fixture rings and raising or lowering the pneumatic arm assembly. When the correct distance is achieved, tighten the fixture rings. Raise the pneumatic arm resting ring until it rests firmly against the bottom fixture ring and tighten it securely.

Mounting Instructions

E
N
G
L
I
S
H

Securing the Bottom of the Column Assembly

3.8

Using a level measuring device, check to see if the round column piece is vertically level. Move the bottom of the column assembly as needed to adjust the level. Once this is accomplished, connect the two remaining (bottom) adjustable arm assemblies to the electrical and pneumatic connector plate on the square column piece. Connect the other end of the arms to the legs of the magnet using the round clamps (figure 3.1).

If you are mounting a sample changer on a 600 MHz magnet that has **Vibration Dampers**, then the bottom of the column sets on the support as pictured in figure 3.1. To position the column so that it is vertically level, turn the adjustable arms to move the column. When the support is used, the bottom adjustable arms are not required.

When the column is fully assembled, place the warning triangle sticker shown at the left, on the upper front half of the column assembly. This warning indicates

that hands and objects should be kept out of the path of the pneumatic arm (refer to chapter 2, section 4).

- í **Before using the sample changer in automatic mode, perform the fine adjustment procedure for the pincher as described in the B-ACS Installation Manual.**


Connecting the Hoses and Cables

3.9

Slip the pneumatic hoses and electrical cables from the arm to the cabinet through the flexible rubber sleeve down, through the column to the side opening into the cabinet (figure 3.4). Connect the pneumatic hoses and the electrical cables from the pneumatic arm assembly, and from the connector panel on the bottom of the column assembly (figure 3.4) to the pneumatic hose connection assembly and electrical sensor connector assembly in the cabinet (figure 3.7, 3.8, 3.9).

Important: All of the electrical cables, electrical connections, air hoses and air hose connections are number-coded. Plug the numbered cable/hose into the corresponding connection (for example: cable # 6 to connection # 6 on the Sensor Connector Assembly).

Figure 3.8. Rear View of the Automatic Sample Changer

Connect the 230V main power cable from an electrical outlet to the Cable and Pneumatic Connection Panel (figure 3.10). Connect the Connection Panel to Cabinet end of the 230V power cable to the Mains Connection as shown in figure 3.8.


Connect cable # 5 to LISH 5 on the Cable and Pneumatic Hose Connection Panel (figure 3.4 and 3.10). Connect the other end of the cable to the Light Barrier for the Shim System as shown in figure 3.2. For magnets with the BOSS 1 and BOSS

Mounting Instructions

2 Shim Systems the cable is connected directly to the BSMS SLCB board, Sample Control connector.


Connect the RS232 connector and cable from the computer to the Remote Control 25-pin female connection on the Cable and Pneumatic Hose Connection Panel (figure 3.10).

Figure 3.9. Partial Front View of Column and Cabinet

E
N
G
L
I
S
H

Figure 3.10. Cable and Pneumatic Hose Connection Panel

Description of Input and Output Locations

3.10

Refer to figure 3.10.

Input: Mains connection - The 230V Connector

Air Input - Air Inlet

Light Barrier Shim System Input

Output: Sample lift air output - Sample Lift. Parallel with buffer.

Input and Output: RS232 Cable connection to computer - Remote Control

Mounting Instructions

E
N
G
L
I
S
H

Special Tools

4

Special Tools for the Automatic Sample Changer

4.1

Table 4.1. Tools for mounting the Automatic Sample Changer

Part Number	Part Name
H1537	Print Board Puller
10499	Open End Wrench 12/13
14667	Open End Wrench 13/17
2920	Open End Wrench 14/17
2850	Allen Key Wrench 2.5 mm CU-BE
10500	Allen Key Wrench 3.0 mm CU-BE
10501	Allen Key Wrench 4.0 mm CU-BE
10502	Allen Key Wrench 5.0 mm CU-BE
14456	Allen Key Wrench 6.0 mm CU-BE
2939	Allen Key Wrench with handle 5 mm
H5042	B-ACS Grease-Box for Vertical Assembly (light- brown color)
H5043	B-ACS Grease-Box for Horizontal Assembly (clear-white color)

General

5.1

During the operation of the sample changer, two kinds of error messages can occur:

- Errors requiring user intervention
- Errors in the communication between the sample changer and the host computer

These two kinds of failures will be treated differently, and are explained in detail in the following section.

Error Messages

5.2

Generally, an error message consists of two parts:

- A two-digit error code
- An explicit error message

A RUB character (decimal 255) will precede any error message transmitted with the SIO channel.

Errors Requiring User Intervention

5.3

1. Insufficient Air Pressure

The actual air pressure is below the allowed limit of 3 bar.

2. Downward Motion Failed

The pincher did not reach its lower vertical position within 10 seconds.

3. Upward Motion Failed

The pincher did not reach its upper vertical position within 10 seconds.

4. Outward Motion Failed

The pincher did not reach its outer position at the magnet within 30 seconds.

5. Inward Motion Failed

The pincher did not reach its inner magazine position within 30 seconds.

6. Pincher Opening Failed

The pincher did not open.

7. Pincher Closing Failed

The pincher did not close.

8. Carrousel Motion Failed

The magazine did not move.

9. Carrousel Position Undefined

The magazine does not lock in the changing position.

10. Sample Holder Not Empty

The magazine position in which the sample should be inserted is already occupied.

11. Sample Detect at Magnet Failed

The optical detector at the magnet sees no sample held by the pincher.

12. Sample Detect at Carrousel Failed

The optical detector at the magazine sees no sample.

13. Sample Detect at Magnet Failed

The optical detector at the magnet sees no sample within 30 seconds of lift on.

14. Sample Grasping Failed

The sample could not be removed.

16. CDW Detection Failed

The number read from the internal code wheel is not within the allowed range between 1 and 60.

21. CDW MISADJUSTED FOR POSITION No: XX

Code XX could not be read from the code wheel.

22. No Free Magazine Position

No free magazine position could be found for the sample in the pincher.

23. Sample Missing

No sample could be found in the present magazine position.

25. Failure of Spinning Device

The mechanism of the bar code reader is no in its required idle position.

26. Arm Positioning Failed

The exact outer magazine position (61 to 120) could not be reached.

If one of the above mentioned failures is observed the program branches off into an error correction routine. It can be restarted by pressing the CONT button or by giving the command „CO” via the RS 232C link. As an alternative the command „HO” can be used to move the sample changer into its „HOME” position.

Failures in the communication between the Sample Changer and the Host Computer can occur if an incorrect command is given by the host computer, or if unconnected hardware is addressed. Following such an error message, the sample changer will accept a new command and not enter the error mode. Error messages are not displayed but are sent back to the host computer. If any of these errors occur consult your operating instructions, or a BRUKER service representative.

50. Bar Code Reader Not Present

An attempt was made to call a command related to the bar code reader even though this is not installed.

51. Invalid Command

The received command was invalid.

52. Invalid Parameter

The command contained an invalid or incorrect parameter.

53. Timeout RS232 ASPECT 3000

A time-out error of the RS232 host TTY has occurred.

54. Horizontal Optic Not Present

The test program for the horizontal optic was called, even though no option is installed on the B-ACS 60.

Operator Maintenance

6

Greasing the Horizontal Arm

6.1

The horizontal arm should be lubricated bi-annually (at least) with B-ACS grease P/N H5043 (has a clear-white color) to allow for smooth operation.

Greasing the Vertical Cylinder

6.2

The moving parts of the vertical cylinder should be lubricated bi-annually (at least) with B-ACS grease P/N H5042 (has a light brown color) to allow for smooth up and down movement.

Arm Adjustment

6.3

The arm should be adjusted initially, and rechecked periodically to ensure that the samples are being properly grasped by the pincher. Refer to the B-ACS Installation Manual, Settings and Adjustments.

Cabinet Adjustment

6.4

The cabinet will need adjusting if the sample, when approaching the magazine belt, is not centered perfectly in the sample holder of the magazine belt. If the inward position of the vertical cylinder is not correct, then refer to the B-ACS Installation Manual, Settings and Adjustments. If the spinner tip edge gets too close to the sample holder (less than 5 mm) move the cabinet downwards on the column (refer to the B-ACS Installation Manual, Settings and Adjustments). The correct distance between the spinner top edge and the magazine belt should be approximately 5 mm.

Vertical Cylinder

6.5

The vertical position of the vertical cylinder is factory set and does not need adjusting. The horizontal position can be adjusted by loosening the two screws under the bottom support bracket of the vertical cylinder (refer to the B-ACS Installation Manual, Settings and Adjustments).

Magazine Belt Tension Adjustment

6.6

The adjustment mechanism under the right carousel wheel serves for setting the tension of the magazine belt. The tension should be increased only when the belt filled with samples hangs more than 5 to 8 mm under the wheel line on the strip between the two wheels. Excessive tension will damage belt and wheel bearings.

Pneumatic Cylinder Setting

6.7

The speed of both vertical and horizontal cylinder can be varied by regulating the outlet air flows of the cylinder supply connections. Refer to the B-ACS Installation Manual, Settings and Adjustments for details.

Cleaning the Inside of the Vertical Cylinder

6.8

When an „Upward Motion Failed“ or „Downward Motion Failed“ error message occurs, it may be necessary to clean the inside of the Vertical Cylinder. This is accomplished by: turning the sample changer off, disconnecting the sensor and hydraulic lines to the vertical cylinder, marking the current position of the cylinder with i.e. a pencil, removing the four screws and cover plate located just above the pincher, and sliding the cylinder out. To clean, wipe the inside of the cylinder with a clean cloth and alcohol. Apply vertical grease P/N H5042 (light brown color) to the cylinder and replace the cylinder, cover and screws and return the cylinder to its original position (that you marked earlier) before tightening. Once you have tightened the screws plug in the sensor and hydraulic lines (for the vertical cylinder) and turn on the sample changer

Technical Data

7

Equipment Identification

7.1

Equipment type: B-ACS

Part Number: B-ACS 60 Sample Changer: H1080
 B-ACS60/600 Sample Changer OXF: H11080
 B-ACS 120 Sample Changer: H800
 B-ACS 60/600 Sample Changer SAG: H5895
 B-ACS 120/600 Sample Changer: H1800

Power Supply Requirements

7.2

110/230V ~, 50/60 Hz
Current carrying capacity: 0.2A

Fuse Protection

7.3

2 x 0.4 AT (slow-blow fuse)

Air Requirements

7.4

3.5-4.5 bar (50-64 psi)
The air source must be **clean, dry and oil-free**.

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