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```
$ ls -l /bin/sh
lrwxrwxrwx 1 root root 4 Mar  1  2012 /bin/sh -> dash
$ sudo dpkg-reconfigure dash
[sudo] password for sermnadmin:
Removing 'diversion of /bin/sh to /bin/sh.distrib by dash'
Adding 'diversion of /bin/sh to /bin/sh.distrib by bash'
Removing 'diversion of /usr/share/man/man1/sh.1.gz to
/usr/share/man/man1/sh.distrib.1.gz by dash'
Adding 'diversion of /usr/share/man/man1/sh.1.gz to
/usr/share/man/man1/sh.distrib.1.gz by bash'
$ ls -l /bin/sh
lrwxrwxrwx 1 root root 4 Aug 24 11:20 /bin/sh -> bash
```

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```
Cannot read FCU BBIS data: cannot open
/opt/topspin1.3/conf/instr/spect/bbis_fcu:
error in "iiconf"
```

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- HPPR Preamplifier 1: `tty10`
- ACB Amplifier Control Board: `tty04`
- RX22 Receiver: `tty10`
- BSMS Smart Magnetc Control System: `tty02`

- Lock Signal: tty03
- LC-NMR Software HyStar: tty07
- VTU Variable Temperature Unit: tty05

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```
Warning: The 2H amplifier in the BSM uses router output 3 which is selected
by the 0.5W 1H Amplifier
```

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```
Operation not permitted
Cannot make link
link must be done manually with
ln /opt/topspin1.3/prog/mod/go4 /opt/topspin1.3/prog/mod/acq_control
```

i

```
Operation not permitted
Cannot make link
link must be done manually with
ln /opt/topspin1.3/prog/wobble/pp_drx_X
/opt/topspin1.3/prog/wobble/pulsprog_X
```

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rpm package file by Bruker	lib files included	Debian package
bruker-compatlibs-1.3-6.el6.brk.x86_64.rpm	libcom_err.so	✔ comerr-dev
—	libcrypto.so	✔ libssl-dev
—	libssl.so	✔ libssl-dev

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```
# Monitor Bruker Topspin spectrometer configuration folder to
# find out which process deletes the bbis_fcw and bbis_rcu files
#
# -w watch /opt/topspin1.3/conf/instr
# -p warx watch for write, attribute change, execute, read events
-w /opt/topspin1.3/conf/instr/ -a exit,always -p warx -k TOPSPIN-CONF

# Monitor the bbis_fcw and bbis_rcu Bruker Topspin spectrometer
# configuration files to find out which process deletes them
#
# -w watch /opt/topspin1.3/conf/instr/spect/bbis_fcw
# -w watch /opt/topspin1.3/conf/instr/spect/bbis_rcu
# -p warx watch for write, attribute change, execute, read events
-w /opt/topspin1.3/conf/instr/spect/bbis_fcw -a exit,always -p warx -k
TOPSPIN-BBIS-FCU
-w /opt/topspin1.3/conf/instr/spect/bbis_rcu1 -a exit,always -p warx -k
TOPSPIN-BBIS-RCU

# Monitor delete files & folders by 64-bit and 32-bit applications
#
# -F arch=b32 or -F arch=b64 monitor 32- or 64-bit applications
# -S unlink watch for unlink (aka delete) system events
```

```
# -S rmdir watch for remove directory system events
# -k TOPSPIN-CONF

-a exit,always -F arch=b64 -S unlink -S unlinkat -S rmdir -k TOPSPIN-DEL-64
-a exit,always -F arch=b64 -S unlink -S unlinkat -S rmdir -k TOPSPIN-DEL-64
```

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```
09:54:52.963 19230 error msg: hconfserver
Operation not permitted
Cannot make link
link must be done manually with
ln /opt/topspin1.3/prog/mod/go4
/opt/topspin1.3/prog/mod/go
09:55:07.362 19230 error msg: hconfserver
Operation not permitted
Cannot make link
link must be done manually with
ln /opt/topspin1.3/prog/mod/go4
/opt/topspin1.3/prog/mod/acq_control
09:55:11.440 19230 error msg: hconfserver
Operation not permitted
Cannot make link
```

```

link must be done manually with
ln /opt/topspin1.3/prog/wobble/pp_drx_X
/opt/topspin1.3/prog/wobble/pulsprog_X

```

Per resoldre'l creo un bash-script `fix_topspin13.sh`

```

#!/bin/sh
#*****
#
# Copyright (c) 2015
# SeRMN, Universitat Autònoma de Barcelona
#
#*****
#           DESCRIPTION
#
# When TopSpin 1.3 runs on a Debian GNU/Linux, the 'cf' command
# fails to create several hard links needed to run the software,
# and to (re)create the program 'shimcntl'.
#
# To fix this problem,
# 1. run 'cf' and ignore all warnings and errors about not
#    being able to create hardlinks.
# 2. run this script. It requires root permissions to create, hence,
#    this script upon running the 'cf' command.
#
#*****

echo ""

# link pulsprog_X to pp_drx_X
echo "creating file '/opt/topspin1.3/prog/wobble/pulsprog_X'"
ln -fv /opt/topspin1.3/prog/wobble/pp_drx_X
/opt/topspin1.3/prog/wobble/pulsprog_X
ls -li /opt/topspin1.3/prog/wobble/pp_drx_X
/opt/topspin1.3/prog/wobble/pulsprog_X
echo ""

# link go to go4
echo "creating file '/opt/topspin1.3/prog/mod/go'"
ln -fv /opt/topspin1.3/prog/mod/go4 /opt/topspin1.3/prog/mod/go
ls -li /opt/topspin1.3/prog/mod/go4 /opt/topspin1.3/prog/mod/go
echo ""

# link acq_control to go4
echo "creating file '/opt/topspin1.3/prog/mod/acq_control'"
ln -fv /opt/topspin1.3/prog/mod/go4 /opt/topspin1.3/prog/mod/acq_control
ls -li /opt/topspin1.3/prog/mod/go4 /opt/topspin1.3/prog/mod/acq_control
echo ""

```

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RCU1: DRAM = 4096 kByte, SRAM = 256 kByte, FIFO = 4 kByte
Digitizers and Filters connected to RCU1:
- HADC/2 ADC937

Gcu: GCU1

Router: 1 Avance-Router

Amplifier Control Board: connected to spect:/dev/tty04

- Warning: The 2H amplifier in the BSMS uses router output 3 which is selected by the 0.5W 1H Amplifier

Amplifiers (ACB configuration):

R0	Amplifier	Module	Nucleus	Power/W	Switchbox	PN/SN/FW
1	1	1	X	300	Y vCB=de	W1301865/564/20000509
2	1	2	1H	50	Y --""--	W1301865/564/20000509
3	2	1	2H	20	Y iCB=00	-/-/-
5	3	1	X	300	vCB=ff	W1301840/1340/19990407

ASU's (BBIS configuration):

AQ-Rack position	MOD	MULT	connected to
6	2	1	FCU-1
6	2	1	FCU-2
5	2	1	FCU-3

(Note: AQ-Rack positions are counted from right)

BSMS: device connected to spect:/dev/tty02

- absolute lockshift frequency = 13227000 Hz
- 2H Amplifier: connected to router output 3

BPSU: device connected to spect:/dev/tty07

- use BACS air = no
- BACS sx delay = 10 s
- Sample Rail fast changer mode = no

Synthesizers :

- PTS-620DL 1MHZ : for F1/F2, 620 MHz max, doubling frequency = 309.00MHz
- PTS-620SL 1MHZ : for F3, 620 MHz max, doubling frequency = 309.00MHz

Lock: device connected to spect:/dev/tty03

Eurotherm VTU: device connected to spect:/dev/tty05

Preamplifiers :

HPPR: - HPPR/2 preamplifier connected to spect:/dev/tty10

Type : HPPR/2

Controller: Application firmware = A0

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```
# cd /home/sermnuab/  
# mkdir data  
# cd data  
# rsync -Hvax --progress --stats /home/sermnadmin/Downloads/Topspin-1.3@old-  
computer/opt_topspin_data/ .  
[...]  
Number of files: 276808  
Number of files transferred: 240678  
Total file size: 12626980788 bytes  
Total transferred file size: 12626980788 bytes  
Literal data: 12626980788 bytes  
Matched data: 0 bytes  
File list size: 4490020  
File list generation time: 0.001 seconds  
File list transfer time: 0.000 seconds  
Total bytes sent: 12643300746  
Total bytes received: 4717413
```

```
sent 12643300746 bytes received 4717413 bytes 10677938.50 bytes/sec
total size is 12626980788 speedup is 1.00
```

```
# chown root:root ../data
# chown -R sermnuab:sermnuab
# chmod 755 ../data
# ls -ld ../data
drwxr-xr-x 28 root root 4096 Sep 22 16:33 ../data
```

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```
No acquisition hardware available for: 192.168.2.7
```

(Original message = `NO_ACQUISITION_AVAILABLE1` 192.168.2.7)

```
=====
Wed Sep 23 16:30:56 CEST 2015
Topspin Version 1.3 (of July 13 2009)
```

```
No acquisition running.
ConfigProp.get: ACQ_DISCFUNC
de.bruker.nmr.jutil.except.MfrException: No acquisition running.
ConfigProp.get: ACQ_DISCFUNC
  at de.bruker.nmr.prsc.cpr.ConfigProp.get(ConfigProp.java:91)
  at de.bruker.nmr.sc.acq.AcqUtil.getAcqFileName(AcqUtil.java:53)
  at de.bruker.nmr.sc.acq.Acq.checkAcquisitionRunning(Acq.java:591)
  at de.bruker.nmr.sc.acq.Acq.showAcqWindow(Acq.java:631)
  at sun.reflect.NativeMethodAccessorImpl.invoke0(Native Method)
  at
sun.reflect.NativeMethodAccessorImpl.invoke(NativeMethodAccessorImpl.java:39
)
  at
sun.reflect.DelegatingMethodAccessorImpl.invoke(DelegatingMethodAccessorImpl
.java:25)
  at java.lang.reflect.Method.invoke(Method.java:585)
  at
de.bruker.nmr.mfw.base.CmdThreadImpl.doJmethod(CmdThreadImpl.java:535)
  at de.bruker.nmr.mfw.base.CmdThreadImpl.exec(CmdThreadImpl.java:439)
  at de.bruker.nmr.mfw.base.CmdThreadImpl.run(CmdThreadImpl.java:393)
```

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```
$ for PROG in /opt/topspin1.3/prog/mod/*;
do
    echo $PROG;
    /opt/topspin1.3/topspin -e ldd $PROG | grep "not found";
done
```

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```
/opt/topspin1.3/prog/mod/xwinplot.mod
libdpstk.so.1 => not found
libdps.so.1 => not found
```

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- https://www.bruker.com/cgi-bin/bkb/show_bug.cgi?id=10711
- https://www.bruker.com/cgi-bin/bkb/show_bug.cgi?id=8259
- https://www.bruker.com/cgi-bin/bkb/show_bug.cgi?id=5956
- https://www.bruker.com/cgi-bin/bkb/show_bug.cgi?id=2375

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