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```
$ ls -l /bin/sh
lrwxrwxrwx 1 root root 4 Mar  1 2012 /bin/sh -> dash
$ sudo dpkg-reconfigure dash
[sudo] password for sermnadmin:
Removing 'diversion of /bin/sh to /bin/sh.distrib by dash'
Adding 'diversion of /bin/sh to /bin/sh.distrib by bash'
Removing 'diversion of /usr/share/man/man1/sh.1.gz to
/usr/share/man/man1/sh.distrib.1.gz by dash'
Adding 'diversion of /usr/share/man/man1/sh.1.gz to
/usr/share/man/man1/sh.distrib.1.gz by bash'
$ ls -l /bin/sh
lrwxrwxrwx 1 root root 4 Aug 24 11:20 /bin/sh -> bash
```

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```
Cannot read FCU BBIS data: cannot open
/opt/topspin1.3/conf/instr/spect/bbis_fcu:
error in "iiconf"
```

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- HPPR Preamplifier 1: `tty10`
- ACB Amplifier Control Board: `tty04`
- RX22 Receiver: `tty10`

- BSMS Smart Magnetic Control System: tty02
- Lock Signal: tty03
- LC-NMR Software HyStar: tty07
- VTU Variable Temperature Unit: tty05

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```
Warning: The 2H amplifier in the BMSM uses router output 3 which is selected
by the 0.5W 1H Amplifier
```

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```
Operation not permitted
Cannot make link
link must be done manually with
ln /opt/topspin1.3/prog/mod/go4 /opt/topspin1.3/prog/mod/acq_control
```

i

```
Operation not permitted
Cannot make link
link must be done manually with
ln /opt/topspin1.3/prog/wobble/pp_drx_X
/opt/topspin1.3/prog/wobble/pulsprog_X
```

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(ignoring RS device entries)

```
#-----  
# FCU configuration:  
#-----  
board # 1 :  
    MEMORY 0x3f000  (=63 k Words)  
    PAL-Version 0x3 = new version (EC level 07 or more)  
board # 2 :  
    MEMORY 0x3f000  (=63 k Words)  
    PAL-Version 0x3 = new version (EC level 07 or more)  
board # 3 :  
    MEMORY 0x3f000  (=63 k Words)  
    PAL-Version 0x3 = new version (EC level 07 or more)
```

RCU1: DRAM = 4096 kByte, SRAM = 256 kByte, FIFO = 4 kByte

Digitizers and Filters connected to RCU1:

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| rpm package file by Bruker                 | lib files included | Debian package |
|--------------------------------------------|--------------------|----------------|
| bruker-compatlibs-1.3-6.el6.brk.x86_64.rpm | libcom_err.so      | ✔ comerr-dev   |
| —                                          | libcrypto.so       | ✔ libssl-dev   |
| —                                          | libssl.so          | ✔ libssl-dev   |

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```
# Monitor Bruker Topspin spectrometer configuration folder to
# find out which process deletes the bbis_fcu and bbis_rcu files
#
# -w watch /opt/topspin1.3/conf/instr
# -p warx watch for write, attribute change, execute, read events

-w /opt/topspin1.3/conf/instr/ -a exit,always -p warx -k TOPSPIN-CONF

# Monitor the bbis_fcu and bbis_rcu Bruker Topspin spectrometer
# configuration files to find out which process deletes them
#
# -w watch /opt/topspin1.3/conf/instr/spect/bbis_fcu
# -w watch /opt/topspin1.3/conf/instr/spect/bbis_rcu
# -p warx watch for write, attribute change, execute, read events

-w /opt/topspin1.3/conf/instr/spect/bbis_fcu -a exit,always -p warx -k
TOPSPIN-BBIS-FCU
-w /opt/topspin1.3/conf/instr/spect/bbis_rcu1 -a exit,always -p warx -k
TOPSPIN-BBIS-RCU
```

```
# Monitor delete files & folders by 64-bit and 32-bit applications
#
# -F arch=b32 or -F arch=b64 monitor 32- or 64-bit applications
# -S unlink watch for unlink (aka delete) system events
# -S rmdir watch for remove directory system events
# -k TOPSPIN-CONF

-a exit,always -F arch=b64 -S unlink -S unlinkat -S rmdir -k TOPSPIN-DEL-64
-a exit,always -F arch=b32 -S unlink -S unlinkat -S rmdir -k TOPSPIN-DEL-32
```

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```
09:54:52.963 19230 error msg: hconfserver
Operation not permitted
Cannot make link
link must be done manually with
ln /opt/topspin1.3/prog/mod/go4
/opt/topspin1.3/prog/mod/go
09:55:07.362 19230 error msg: hconfserver
Operation not permitted
Cannot make link
link must be done manually with
```

```

ln /opt/topspin1.3/prog/mod/go4
/opt/topspin1.3/prog/mod/acq_control
09:55:11.440 19230 error msg: hconfserver
Operation not permitted
Cannot make link
link must be done manually with
ln /opt/topspin1.3/prog/wobble/pp_drx_X
/opt/topspin1.3/prog/wobble/pulsprog_X

```

Per resoldre'l creo un bash-script *fix\_topspin13.sh*

```

#!/bin/sh
#*****
#
# Copyright (c) 2015
# SeRMN, Universitat Autònoma de Barcelona
#
#*****
# DESCRIPTION
#
# When TopSpin 1.3 runs on a Debian GNU/Linux, the 'cf' command
# fails to create several hard links needed to run the software,
# and to (re)create the program 'shimcntl'.
#
# To fix this problem,
# 1. run 'cf' and ignore all warnings and errors about not
#    being able to create hardlinks.
# 2. run this script. It requires root permissions to create, hence,
#    this script upon running the 'cf' command.
#
#*****

echo ""

# link pulsprog_X to pp_drx_X
echo "creating file '/opt/topspin1.3/prog/wobble/pulsprog_X'"
ln -fv /opt/topspin1.3/prog/wobble/pp_drx_X
/opt/topspin1.3/prog/wobble/pulsprog_X
ls -li /opt/topspin1.3/prog/wobble/pp_drx_X
/opt/topspin1.3/prog/wobble/pulsprog_X
echo ""

# link go to go4
echo "creating file '/opt/topspin1.3/prog/mod/go'"
ln -fv /opt/topspin1.3/prog/mod/go4 /opt/topspin1.3/prog/mod/go
ls -li /opt/topspin1.3/prog/mod/go4 /opt/topspin1.3/prog/mod/go
echo ""

# link acq_control to go4
echo "creating file '/opt/topspin1.3/prog/mod/acq_control'"
ln -fv /opt/topspin1.3/prog/mod/go4 /opt/topspin1.3/prog/mod/acq_control

```

```
ls -li /opt/topspin1.3/prog/mod/go4 /opt/topspin1.3/prog/mod/acq_control  
echo ""
```

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```
RCU1: DRAM = 4096 kByte, SRAM = 256 kByte, FIFO = 4 kByte  
Digitizers and Filters connected to RCU1:  
- HADC/2 ADC937
```

```
Gcu: GCU1
```

```
Router: 1 Avance-Router
```

```
Amplifier Control Board:  connected to spect:/dev/tty04
```

```
- Warning: The 2H amplifier in the BSMS uses router output 3 which is  
          selected by the 0.5W 1H Amplifier
```

```
Amplifiers (ACB configuration):
```

| R0 | Amplifier | Module | Nucleus | Power/W | Switchbox | PN/SN/FW               |
|----|-----------|--------|---------|---------|-----------|------------------------|
| 1  | 1         | 1      | X       | 300     | Y vCB=de  | W1301865/564/20000509  |
| 2  | 1         | 2      | 1H      | 50      | Y --" "-- | W1301865/564/20000509  |
| 3  | 2         | 1      | 2H      | 20      | Y iCB=00  | -/-/-                  |
| 5  | 3         | 1      | X       | 300     | vCB=ff    | W1301840/1340/19990407 |

ASU's (BBIS configuration):

| AQ-Rack position | MOD | MULT | connected to |
|------------------|-----|------|--------------|
| 6                | 2   | 1    | FCU-1        |
| 6                | 2   | 1    | FCU-2        |
| 5                | 2   | 1    | FCU-3        |

(Note: AQ-Rack positions are counted from right)

BSMS: device connected to spect:/dev/tty02

- absolute lockshift frequency = 13227000 Hz
- 2H Amplifier: connected to router output 3

BPSU: device connected to spect:/dev/tty07

- use BACS air = no
- BACS sx delay = 10 s
- Sample Rail fast changer mode = no

Synthesizers :

- PTS-620DL 1MHZ : for F1/F2, 620 MHz max, doubling frequency = 309.00MHz
- PTS-620SL 1MHZ : for F3, 620 MHz max, doubling frequency = 309.00MHz

Lock: device connected to spect:/dev/tty03

Eurotherm VTU: device connected to spect:/dev/tty05

Preamplifiers :

HPPR: - HPPR/2 preamplifier connected to spect:/dev/tty10

Type : HPPR/2

Controller: Application firmware = A0

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```

Total file size: 12626980788 bytes
Total transferred file size: 12626980788 bytes
Literal data: 12626980788 bytes
Matched data: 0 bytes
File list size: 4490020
File list generation time: 0.001 seconds
File list transfer time: 0.000 seconds
Total bytes sent: 12643300746
Total bytes received: 4717413

sent 12643300746 bytes  received 4717413 bytes  10677938.50 bytes/sec
total size is 12626980788  speedup is 1.00

# chown root:root ../data
# chown -R sermnuab:sermnuab
# chmod 755 ../data
# ls -ld ../data
drwxr-xr-x 28 root root 4096 Sep 22 16:33 ../data

```

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```
# diff topspin topspin.Orig-20150824
62c62
<     export
LD_LIBRARY_PATH="$szProgPath"/pvshlib:"$szProgPath"/shlib:$LD_LIBRARY_PATH
---
>     export LD_LIBRARY_PATH="$szProgPath"/pvshlib:"$szProgPath"/shlib
```

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```
# MCE - 20150925 - Library path environment variables added to allow
# compiling TopSpin 1.3 AU programs in newer GNU/Linux versions

export PATH=$PATH:/opt/topspin1.3/gnu/bin
export LIBRARY_PATH=/usr/lib32:$LIBRARY_PATH
export LD_LIBRARY_PATH=/usr/lib32:/lib32:$LD_LIBRARY_PATH
```

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```
/opt/topspin1.3/exp/stan/nmr/au/src/coiltemp:42:5: error: conflicting types
for 'coiltemp_program'
/opt/topspin1.3/exp/stan/nmr/au/src/coiltemp:43:1: note: an argument type
that has a default promotion can't match an empty parameter name list
declaration
/opt/topspin1.3/exp/stan/nmr/au/src/coiltemp:31:9: note: previous implicit
declaration of 'coiltemp_program' was here
```

```
exit 1
/opt/topspin1.3/exp/stan/nmr/au/src/decra: In function 'AU_program':
/opt/topspin1.3/exp/stan/nmr/au/src/decra:112:12: error: invalid storage
class for function 'DECRA_program'
/opt/topspin1.3/exp/stan/nmr/au/src/decra: At top level:
/opt/topspin1.3/exp/stan/nmr/au/src/decra:293:12: error: static declaration
of 'DECRA_program' follows non-static declaration
/opt/topspin1.3/exp/stan/nmr/au/src/decra:115:9: note: previous implicit
declaration of 'DECRA_program' was here
exit 1
/opt/topspin1.3/exp/stan/nmr/au/src/disconv: In function 'AU_program':
/opt/topspin1.3/exp/stan/nmr/au/src/disconv:128:12: error: invalid storage
class for function 'DISNMR_program'
/opt/topspin1.3/exp/stan/nmr/au/src/disconv: At top level:
/opt/topspin1.3/exp/stan/nmr/au/src/disconv:878:12: error: static
declaration of 'DISNMR_program' follows non-static declaration
/opt/topspin1.3/exp/stan/nmr/au/src/disconv:131:9: note: previous implicit
declaration of 'DISNMR_program' was here
exit 1
/opt/topspin1.3/exp/stan/nmr/au/src/disinfo: In function 'AU_program':
/opt/topspin1.3/exp/stan/nmr/au/src/disinfo:106:12: error: invalid storage
class for function 'DISNMR_program'
/opt/topspin1.3/exp/stan/nmr/au/src/disinfo: At top level:
/opt/topspin1.3/exp/stan/nmr/au/src/disinfo:832:12: error: static
declaration of 'DISNMR_program' follows non-static declaration
/opt/topspin1.3/exp/stan/nmr/au/src/disinfo:109:9: note: previous implicit
declaration of 'DISNMR_program' was here
exit 1
/opt/topspin1.3/exp/stan/nmr/au/src/lcsino:46:7: error: conflicting types
for 'lcsino'
/opt/topspin1.3/exp/stan/nmr/au/src/lcsino:40:7: note: previous implicit
declaration of 'lcsino' was here
exit 1
/opt/topspin1.3/exp/stan/nmr/au/src/multi_zgvd: In function 'AU_program':
/opt/topspin1.3/exp/stan/nmr/au/src/multi_zgvd:47:13: error: invalid storage
class for function 'PrintExpTime'
exit 1
/opt/topspin1.3/exp/stan/nmr/au/src/multi_zgvt: In function 'AU_program':
/opt/topspin1.3/exp/stan/nmr/au/src/multi_zgvt:54:13: error: invalid storage
class for function 'PrintExpTime'
exit 1
/opt/topspin1.3/exp/stan/nmr/au/src/multicmd: In function 'AU_program':
/opt/topspin1.3/exp/stan/nmr/au/src/multicmd:39:13: error: invalid storage
class for function 'PrintExpTime'
exit 1
/opt/topspin1.3/exp/stan/nmr/au/src/multiexpt: In function 'AU_program':
/opt/topspin1.3/exp/stan/nmr/au/src/multiexpt:29:13: error: invalid storage
class for function 'PrintExpTime'
exit 1
```

```
/opt/topspin1.3/exp/stan/nmr/au/src/multizg: In function 'AU_program':
/opt/topspin1.3/exp/stan/nmr/au/src/multizg:38:13: error: invalid storage
class for function 'PrintExpTime'
exit 1
/opt/topspin1.3/exp/stan/nmr/au/src/noemult: In function 'AU_program':
/opt/topspin1.3/exp/stan/nmr/au/src/noemult:53:13: error: invalid storage
class for function 'PrintExpTime'
exit 1
/opt/topspin1.3/exp/stan/nmr/au/src/pulsecalib: In function 'AU_program':
/opt/topspin1.3/exp/stan/nmr/au/src/pulsecalib:123:16: error: invalid
storage class for function 'pulsecalib'
/opt/topspin1.3/exp/stan/nmr/au/src/pulsecalib: At top level:
/opt/topspin1.3/exp/stan/nmr/au/src/pulsecalib:290:1: error: static
declaration of 'pulsecalib' follows non-static declaration
/opt/topspin1.3/exp/stan/nmr/au/src/pulsecalib:124:12: note: previous
implicit declaration of 'pulsecalib' was here
exit 1
/opt/topspin1.3/exp/stan/nmr/au/src/setproj:149:5: error: conflicting types
for 'sorter'
/opt/topspin1.3/exp/stan/nmr/au/src/setproj:35:5: note: previous declaration
of 'sorter' was here
exit 1
/opt/topspin1.3/exp/stan/nmr/au/src/simplex: In function 'AU_program':
/opt/topspin1.3/exp/stan/nmr/au/src/simplex:262:12: error: invalid storage
class for function 'SIMPLEX_program'
/opt/topspin1.3/exp/stan/nmr/au/src/simplex: At top level:
/opt/topspin1.3/exp/stan/nmr/au/src/simplex:468:12: error: static
declaration of 'SIMPLEX_program' follows non-static declaration
/opt/topspin1.3/exp/stan/nmr/au/src/simplex:265:9: note: previous implicit
declaration of 'SIMPLEX_program' was here
exit 1
```

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```
ILLEGAL TERMINATION of process 28469
due to signal 11
```

potser pels motius explicats [aquí](#): *ldd with 32-bit binaries on 64-bit system - runs but won't report.*  
</note>

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```
$ for PROG in /opt/topspin1.3/prog/mod/*;
do
    echo $PROG;
    /opt/topspin1.3/topspin -e ldd $PROG | grep "not found";
done
```

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```
/opt/topspin1.3/prog/mod/xwinplot.mod
libdpstk.so.1 => not found
libdps.so.1 => not found
```

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- [https://www.bruker.com/cgi-bin/bkb/show\\_bug.cgi?id=10711](https://www.bruker.com/cgi-bin/bkb/show_bug.cgi?id=10711)
- [https://www.bruker.com/cgi-bin/bkb/show\\_bug.cgi?id=8259](https://www.bruker.com/cgi-bin/bkb/show_bug.cgi?id=8259)
- [https://www.bruker.com/cgi-bin/bkb/show\\_bug.cgi?id=5956](https://www.bruker.com/cgi-bin/bkb/show_bug.cgi?id=5956)
- [https://www.bruker.com/cgi-bin/bkb/show\\_bug.cgi?id=2375](https://www.bruker.com/cgi-bin/bkb/show_bug.cgi?id=2375)

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