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## 2017-05-18 - TREBALL EN MARXA / WORK IN PROGRESS

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- [Debian multiarch support](#)
- [Multiarch HOWTO for users](#). User-oriented info on how to use multiarch: configuration, installation, cross-building.

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- Install Bruker NMR Software
- Install GPL licensed GNU GCC 4.9.1 compiler (recommended)

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```
Note: Using init directory /etc/init.d
Note: Using Systemd for service control
```

```
/usr/diskless/prog/bin/setfirewall.sh: install firewall rules for Spectrometer Operating System (show only)
command: iptables-restore < /etc/sysconfig/iptables
saving of iptables not necessary
command: ip6tables-restore < /etc/sysconfig/ip6tables
saving of ip6tables not necessary
/usr/diskless/prog/bin/setfirewall.sh: terminated with status OK
```

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```
*** Error ***

Could not show firewall rules for MICS:
(command: /bin/bash /opt/Bruker/mics/setfirewall.sh show)

/opt/Bruker/mics/service-lib.sh: line 241: rpm: command not found
Fatal error: is not a number
/opt/Bruker/mics/setfirewall.sh: terminated with status 3
```

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```
*** Warning ***

Cannot create user account 'flexlm' and/or group 'flexlm':
Cannot create user account flexlm: useradd flexlm : useradd: group flexlm
exists - if you want to add this user to that group, use -g.
```

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```
installation of DISKLESS 3.5 finished>
>> installnrmr terminated at 2017-05-18 11:51:47 +0200
>> STEP post terminated for mod=hostacq with result = NEXT
>> STEP post started for mod=ipso at 2017-05-18 11:51:47 +0200 ...
>> cd /tftpboot
>> Status: Running post-installation scripts for Ipso0S
>> ::file::write (version.txt, -log)
>> Status: Post processing for Ipso0S ...
>> Execute cmd: LAUNCH=</bin/bash /tftpboot/install-services >
>> Execute cmd: options log=0 canfail=0 cont=0 argc=2>
>> Execute cmd: pwd=</tftpboot>

*** ERROR: ***

Could not enable service tftp: not in xinetd.conf

Failed to get unit file state for dhcpd.service: No such file or directory
Failed to execute operation: No such file or directory
(systemctl enable dhcpd.service) [FAILED]
sed: can't read /usr/lib/systemd/system/dhcpd.service: No such file or
directory
/tftpboot/service-lib.sh: line 1272: restorecon: command not found

*** ERROR: ***

restorecon -v /etc/systemd/system/dhcpd.service

/tftpboot/service-lib.sh: line 1272: restorecon: command not found

*** ERROR: ***

restorecon -v /etc/dhcpd.conf

Failed to start dhcpd.service: Unit dhcpd.service failed to load: Invalid
argument. See system logs and 'systemctl status dhcpd.service' for details.
(systemctl start dhcpd.service) [FAILED]
(systemctl status dhcpd.service) [FAILED]

Output:
Note: Using init directory /etc/init.d
Note: Using Systemd for service control

*** WARNING: ***

network interface eth1 not found
```

```
/tftpboot/install-services: install services: tftp dhcpd
systemctl enable dhcpd.service ... command: echo > /tftpboot/dhcpd.service
command: cp -af /tftpboot/dhcpd.service /etc/systemd/system/dhcpd.service
command: restorecon -v /etc/systemd/system/dhcpd.service
command: cp -af /tftpboot/dhcpd.conf /etc/dhcpd.conf
command: restorecon -v /etc/dhcpd.conf
systemctl -q daemon-reload ... [OK]
systemctl start dhcpd.service ... systemctl status dhcpd.service ... *
dhcpd.service
  Loaded: error (Reason: Invalid argument)
  Active: inactive (dead)

May 18 11:51:47 cie-55-31 systemd[1]: dhcpd.service lacks both ExecStart=
and ExecStop= setting. Refusing.
May 18 11:51:47 cie-55-31 systemd[1]: dhcpd.service lacks both ExecStart=
and ExecStop= setting. Refusing.
/tftpboot/install-services: terminated with error status 1

Saved in: /opt/topspin3.5pl6/install.log
```

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```
during command /tftpboot/install-services
Errors:
/tftpboot/service-lib.sh: line 241: rpm: command not found
*** Fatal error: ***

/tftpboot/install-services: is not a number

Output:
/tftpboot/install-services: terminated with status 99
```

i

```
Error: Execute cmd: exit=99
Error: Execute cmd: error=</usr/diskless/crco_data/service-lib.sh: line 241:
rpm: command not found
*** Fatal error: ***
```

```
/usr/diskless/crco_data/crco-install: is not a number
>
Error: Execute cmd: exit=99
Error: Execute cmd: error=</tftpboot/service-lib.sh: line 241: rpm: command
not found
*** Fatal error: ***

/tftpboot/install-services: is not a number
```

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```
during command installation script (post)
Errors:
./Bruker/service-lib.sh: line 241: rpm: command not found
is not a number
/home/sermnadmin/Downloads/Topspin-3.2.7/Topspin3.2PL7/xwinstall.d/modscript
s/flexlm.post: terminated with status 99

Output:
```

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```
*** Error ***
```

```
*
Could not show firewall rules for Diskless:
(command: /bin/bash /usr/diskless/prog/bin/setfirewall.sh show)

/usr/diskless/prog/bin/service-lib.sh: line 241: rpm: command not found
Fatal error: is not a number
/usr/diskless/prog/bin/setfirewall.sh: terminated with status 3
```

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```
sermadmin@cie-55-31:~$ topspin
TopSpin 3.5 pl 6 - Copyright (C) 2016 Bruker BioSpin GmbH
Installation directory: "/opt/topspin3.5pl6"
Waiting for FLEXlm license

Cannot checkout TopSpin FLEXlm license
Invalid license file syntax
Feature:      TOPSPIN3
License path: /usr/local/flexlm/Bruker/licenses/license.dat
FLEXlm error: -2,413
For further information, refer to the FLEXlm End User Manual,
available at "www.macrovision.com".

The FLEXlm host ID of this machine is 0010182755de

    If your workstation controls the spectrometer,
    you may run the program without license.
    This requires the acquisition server to be ready.

checking connection to spectrometer (may last several minutes)
CPR : 2017-05-19 17:30:51.941 +0200
TopSpin is currently running under
the spectrometer service license.
Please apply to receive a correct
license for your workstation
using the online license form:
https://www.bruker.com/nmr_license_requests.html

error message cprserver: GetSystemPrinterList: dlopen(libcups.so) failed
Exception in thread "main" java.lang.UnsatisfiedLinkError:
/opt/topspin3.5pl6/jre/lib/i386/libawt_xawt.so: libXtst.so.6: cannot open
shared object file: No such file or directory
    at java.lang.ClassLoader$NativeLibrary.load(Native Method)
```

```
at java.lang.ClassLoader.loadLibrary0(ClassLoader.java:1941)
at java.lang.ClassLoader.loadLibrary(ClassLoader.java:1824)
at java.lang.Runtime.load0(Runtime.java:809)
at java.lang.System.load(System.java:1086)
at java.lang.ClassLoader$NativeLibrary.load(Native Method)
at java.lang.ClassLoader.loadLibrary0(ClassLoader.java:1941)
at java.lang.ClassLoader.loadLibrary(ClassLoader.java:1845)
at java.lang.Runtime.loadLibrary0(Runtime.java:870)
at java.lang.System.loadLibrary(System.java:1122)
at java.awt.Toolkit$3.run(Toolkit.java:1636)
at java.awt.Toolkit$3.run(Toolkit.java:1634)
at java.security.AccessController.doPrivileged(Native Method)
at java.awt.Toolkit.loadLibraries(Toolkit.java:1633)
at java.awt.Toolkit.<clinit>(Toolkit.java:1668)
at java.awt.Color.<clinit>(Color.java:275)
at de.bruker.nmr.mfw.base.BSplashScreen.<clinit>(BSplashScreen.java:44)
at de.bruker.nmr.mfw.base.MfwBuilder.programMain(MfwBuilder.java:420)
at bruker.bio.start.topspin.Start.main(Start.java:47)
```

premature Java Virtual Machine termination  
Program is exiting ...

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Ara el TopSpin ja s'executa,

```

sermnadmin@cie-55-31:~$ topspin
TopSpin 3.5 pl 6 - Copyright (C) 2016 Bruker BioSpin GmbH
Installation directory: "/opt/topspin3.5pl6"
Waiting for FLEXlm license

Cannot checkout TopSpin FLEXlm license
Invalid license file syntax

```

```
Feature: TOPSPIN3
License path: /usr/local/flexlm/Bruker/licenses/license.dat
FLEXlm error: -2,413
For further information, refer to the FLEXlm End User Manual,
available at "www.macrovision.com".
```

The FLEXlm host ID of this machine is 0010182755de

If your workstation controls the spectrometer,  
you may run the program without license.  
This requires the acquisition server to be ready.

```
checking connection to spectrometer (may last several minutes)
CPR : 2017-05-19 17:33:57.165 +0200
TopSpin is currently running under
the spectrometer service license.
Please apply to receive a correct
license for your workstation
using the online license form:
https://www.bruker.com/nmr_license_requests.html
```

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IPSO: connected to spectrometer subnet

- TCP/IP address = 149.236.99.90
- IPS0 type: AQS
- Tctrl : 1
- Fctrls: 2
- Gctrl1: with digital preemphasis
- Rctrl : 1
- FREDs : none
- DPP : 1

DRU1: AQS DRU Z100977/00582 ECL 04.00

- TCP/IP address = 149.236.99.89
- Firmware Version = 151014
- DRU1 controls AQS-Rack and HPPR/2

AQS: connected to 149.236.99.89:/dev/tty10

Slot	SBSB	Board				ECL	Name	Description
Number	Addr	Type	HW-VS	FW-VS	ID			
2	0x10	0x72	0x3	AW	R	3.0	REC-1	AQS RXAD400 Z102116/705 ECL 03.00
3	0x34	0x2	0x1		X	5.0	REF-1	REF-400 Reference Board for AQS Receiver
4	0x24	0x11	0x8	AE	S	2.0	SGU-1	AQS SGU/2 400 Z103080/00351 ECL 02.00 FPGA-Vs=AF
5	0x25	0x11	0x8	AE	S	2.0	SGU-2	AQS SGU/2 400 Z103080/00340 ECL 02.00 FPGA-Vs=AF
9	0x36	0x5	0		A	6.0	BLA-1	BLA2BB 150/60 20-400 W1345049/1817 ECL 6
21	0x21	0xcf	0		P	1.0	PSD-1	PSD/3 Power Supply Distribution Board
--	0x20	0xd1	0		B	2.0	MASTER	AQS/3 Chassis
1	0x20	0x7	0		B		MASTER	AQS Rack Master

Logical Receiver/SGU configuration:

Router: 1 AQS-Minirouter

Amplifiers (AQS configuration):

RO	Amplifier	Module	Nucleus	Power/W	Switchbox	Name	PN/SN/FW
1	1	1	B	150	N	BLA2BB	150/60 20-400 W1345049/1817/-
2	1	2	B	60	N	BLA2BB	150/60 20-400 W1345049/1817/-
3	2	1	2H	20	N	2HTX in BSMS	-/-/-

Gradient amplifiers at the spectrometer subnet:

BGA1: BGA\_W1213891\_0040

- TCP/IP address = 149.236.99.93
- Firmware version = 20080206
- Web version = 1.1
- Current limits = 0.0/X, 0.0/Y, -10.0/Z (in A)

BSMS: BSMS/2 connected to ethernet

- TCP/IP address = 149.236.99.20
- ELCB firmware version = 20150423
- SCB channels = 20
- Shim System = BOSS1-SB
- 2HTX Amplifier: available
- Lock: installed

MAS2 control unit: connected to 149.236.99.90:/dev/tty05

- Boot Software : 20031216
- Application Software : 20120329
- Insert Time : 10 s
- Eject Time : 10 s
- Spin Lock Tolerance : 5 Hz

Eurotherm VTU: device connected to 149.236.99.90:/dev/tty04

Preamplifiers :

HPPR2: - HPPR/2 preamplifier connected to 149.236.99.89:/dev/tty10

Type : HPPR/2

Controller: Cover/1 with Application firmware = AQ  
15 LEDs for tuning, 15 LEDs for matching

Module 1 : 1H/2H 1H (virtual 50 Ohm reference: 274.9%/85.9deg)  
PN=Z003959, SN=00418 from 20070319

Module 2 : 1H/2H 2H (virtual 50 Ohm reference: 116.6%/-71.1deg)  
PN=Z003959, SN=00418 from 20070319

Module 3 : XBB19F 2HS  
PN=Z003504, SN=00415 from 20070425

Gradient Controller cable connections

-----  
GCtrl1: connected via DPP to GREAT 1-10 Gradient Amplifier Unit  
W1213891/0040

RF cable connections (assumed)

-----  
SGU1 NORM output -> input 2 of transmitter 1 (BLA2BB 150/60 20-400  
W1345049/1817)

SGU1 AUX output -> input 1 of transmitter 2 (2HTX in BSMS)

SGU2 NORM output -> input 3 of transmitter 1 (BLA2BB 150/60 20-400  
W1345049/1817)

SGU2 AUX output -> TUNE signal input of HPPR

Blanking cable connections (assumed)

-----  
transmitter 1 = BLA2BB 150/60 20-400 W1345049/1817:

- amplifier B-150W uses blanking 1
- amplifier B-60W uses blanking 2

```
transmitter 2 = 2HTX in BSMS:  
- amplifier 2H-20W uses blanking 3
```

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```
atma_ProbeCorba_impl::start exception: Error:can't stat the wobble process  
[NucBB-NucBB- ]
```

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-  `<TOPSPIN_HOME>/exp/stan/nmr/par/users/`
-  `<TOPSPIN_HOME>/exp/stan/nmr/par/teo*` (i `pau*`, `mir*`. Tots es copien a `<TOPSPIN_HOME>/exp/stan/nmr/par/users/`)
-  `<TOPSPIN_HOME>/exp/stan/nmr/lists/pp/users/`
-  `<TOPSPIN_HOME>/exp/stan/nmr/lists/gp/users/` (carpeta buida)
-  `<TOPSPIN_HOME>/exp/stan/nmr/lists/vc/` (carpeta buida)
-  `<TOPSPIN_HOME>/exp/stan/nmr/lists/vd/`
-  `<TOPSPIN_HOME>/exp/stan/nmr/lists/vt/` (carpeta buida)

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```
# 20160422 - MCE - Hide nmrsu user from the login screen
IncludeAll = false
Include = sermnadmin,sermnuab
```

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```
# 20160422 - MCE - Prevent nmrsu user from log-in via ssh
DenyUsers nmrsu
```

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```
root@cie-55-31:/home/sermnadmin/Documents/old-old-hdd-recovery/sdb4-old-old-hdd-opt/topspin/data# rsync -Hvax --progress --stats ./*/opt/topspin3.5pl6/data/[...]
Number of files: 7,573 (reg: 6,405, dir: 1,168)
Number of created files: 7,573 (reg: 6,405, dir: 1,168)
Number of deleted files: 0
Number of regular files transferred: 6,405
Total file size: 459,101,587 bytes
Total transferred file size: 459,101,587 bytes
Literal data: 459,101,587 bytes
Matched data: 0 bytes
File list size: 65,536
File list generation time: 0.001 seconds
File list transfer time: 0.000 seconds
Total bytes sent: 459,608,752
Total bytes received: 128,851

sent 459,608,752 bytes  received 128,851 bytes  54,086,776.82 bytes/sec
total size is 459,101,587  speedup is 1.00
```

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```
root@cie-55-31:/home/sermnadmin/Documents/old-old-hdd-recovery/sdb4-old-old-hdd-opt/sermn-uab# rsync -Hvax --progress --stats ./data /home/sermnuab/
```

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```
root@cie-55-31:/etc# cat /etc/exports
# /etc/exports: the access control list for filesystems which may be
exported
#
#           to NFS clients.  See exports(5).
#
# Example for NFSv2 and NFSv3:
# /srv/homes          hostname1(rw,sync,no_subtree_check)
hostname2(ro,sync,no_subtree_check)
#
# Example for NFSv4:
```

```
# /srv/nfs4          gss/krb5i(rw,sync,fsid=0,crossmnt,no_subtree_check)
# /srv/nfs4/homes   gss/krb5i(rw,sync,no_subtree_check)
#

# SERMN - EXPORTED FOLDERS

# Export data folder to SeRMN data server
/home/sermnuab/data
sermn.uab.es(ro,sync,no_subtree_check,root_squash,all_squash,insecure)

# Export TopSpin 3.5 pl6 lists folder to SeRMN data server
/opt/topspin3.5pl6/exp/stan/nmr/lists
sermn.uab.es(ro,sync,no_subtree_check,root_squash,all_squash,insecure)

# DO NOT MAKE CHANGES BELOW THIS LINE
/usr/diskless/dl_usr          spect(sync,ro)
/usr/diskless/clients/spect   spect(sync,rw,no_root_squash)
```

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```
root@cie-55-31:/etc# exportfs -av
exportfs: /etc/exports [7]: Neither 'subtree_check' or 'no_subtree_check'
specified for export "spect:/usr/diskless/dl
_usr".
  Assuming default behaviour ('no_subtree_check').
  NOTE: this default has changed since nfs-utils version 1.0.x

exportfs: /etc/exports [8]: Neither 'subtree_check' or 'no_subtree_check'
specified for export "spect:/usr/diskless/cl
ients/spect".
  Assuming default behaviour ('no_subtree_check').
  NOTE: this default has changed since nfs-utils version 1.0.x

exporting spect:/usr/diskless/clients/spect
exporting spect:/usr/diskless/dl_usr
exporting sermn.uab.es:/opt/topspin3.5pl6/exp/stan/nmr/lists
exporting sermn.uab.es:/home/sermnuab/data
```

Llavors provo a muntar les carpetes al servidor de dades i no em dona cap error, i es pot accedir al contingut de les carpetes.

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Nuclei	Probe: 1H-BB/19F BBO+ATMA			Console: SGU & Amp
	ATMA / ATMM (Tune & Match)	Power Limit	PROSOL (Pulses)	CORTAB (RF Calibration)
1H	✓	✓	✓	✓
13C	✓	✓	✓	✓
31P	✓	✓	✓	✓
19F	✓	✓	✓	✓
2H	✓	✓	✓	✓

From: <https://sermn.uab.cat/wiki/> - **SeRMN Wiki**

Permanent link: [https://sermn.uab.cat/wiki/doku.php?id=informatica:hp\\_xw4400\\_400sb\\_topspin35&rev=1496408294](https://sermn.uab.cat/wiki/doku.php?id=informatica:hp_xw4400_400sb_topspin35&rev=1496408294)

Last update: **2017/06/02 14:58**