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2017-05-18 - TREBALL EN MARXA / WORK IN PROGRESS

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- [Debian multiarch support](#)
- [Multiarch HOWTO for users](#). User-oriented info on how to use multiarch: configuration, installation, cross-building.

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- Install Bruker NMR Software
- Install GPL licensed GNU GCC 4.9.1 compiler (recommended)

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```
Note: Using init directory /etc/init.d
Note: Using Systemd for service control
```

```
/usr/diskless/prog/bin/setfirewall.sh: install firewall rules for Spectrometer Operating System (show only)
command: iptables-restore < /etc/sysconfig/iptables
saving of iptables not necessary
command: ip6tables-restore < /etc/sysconfig/ip6tables
saving of ip6tables not necessary
/usr/diskless/prog/bin/setfirewall.sh: terminated with status OK
```

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```
*** Error ***

Could not show firewall rules for MICS:
(command: /bin/bash /opt/Bruker/mics/setfirewall.sh show)

/opt/Bruker/mics/service-lib.sh: line 241: rpm: command not found
Fatal error: is not a number
/opt/Bruker/mics/setfirewall.sh: terminated with status 3
```

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```
*** Warning ***

Cannot create user account 'flexlm' and/or group 'flexlm':
Cannot create user account flexlm: useradd flexlm : useradd: group flexlm
exists - if you want to add this user to that group, use -g.
```

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```
installation of DISKLESS 3.5 finished>
>> installnrmr terminated at 2017-05-18 11:51:47 +0200
>> STEP post terminated for mod=hostacq with result = NEXT
>> STEP post started for mod=ipso at 2017-05-18 11:51:47 +0200 ...
>> cd /tftpboot
>> Status: Running post-installation scripts for Ipso0S
>> ::file::write (version.txt, -log)
>> Status: Post processing for Ipso0S ...
>> Execute cmd: LAUNCH=</bin/bash /tftpboot/install-services >
>> Execute cmd: options log=0 canfail=0 cont=0 argc=2>
>> Execute cmd: pwd=</tftpboot>

*** ERROR: ***

Could not enable service tftp: not in xinetd.conf

Failed to get unit file state for dhcpd.service: No such file or directory
Failed to execute operation: No such file or directory
(systemctl enable dhcpd.service) [FAILED]
sed: can't read /usr/lib/systemd/system/dhcpd.service: No such file or
directory
/tftpboot/service-lib.sh: line 1272: restorecon: command not found

*** ERROR: ***

restorecon -v /etc/systemd/system/dhcpd.service

/tftpboot/service-lib.sh: line 1272: restorecon: command not found

*** ERROR: ***

restorecon -v /etc/dhcpd.conf

Failed to start dhcpd.service: Unit dhcpd.service failed to load: Invalid
argument. See system logs and 'systemctl status dhcpd.service' for details.
(systemctl start dhcpd.service) [FAILED]
(systemctl status dhcpd.service) [FAILED]

Output:
Note: Using init directory /etc/init.d
Note: Using Systemd for service control

*** WARNING: ***

network interface eth1 not found
```

```

/tftpboot/install-services: install services: tftp dhcpd
systemctl enable dhcpd.service ... command: echo > /tftpboot/dhcpd.service
command: cp -af /tftpboot/dhcpd.service /etc/systemd/system/dhcpd.service
command: restorecon -v /etc/systemd/system/dhcpd.service
command: cp -af /tftpboot/dhcpd.conf /etc/dhcpd.conf
command: restorecon -v /etc/dhcpd.conf
systemctl -q daemon-reload ... [OK]
systemctl start dhcpd.service ... systemctl status dhcpd.service ... *
dhcpd.service
  Loaded: error (Reason: Invalid argument)
  Active: inactive (dead)

May 18 11:51:47 cie-55-31 systemd[1]: dhcpd.service lacks both ExecStart=
and ExecStop= setting. Refusing.
May 18 11:51:47 cie-55-31 systemd[1]: dhcpd.service lacks both ExecStart=
and ExecStop= setting. Refusing.
/tftpboot/install-services: terminated with error status 1

Saved in: /opt/topspin3.5pl6/install.log

```

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```

during command /tftpboot/install-services
Errors:
/tftpboot/service-lib.sh: line 241: rpm: command not found
*** Fatal error: ***

/tftpboot/install-services: is not a number

Output:
/tftpboot/install-services: terminated with status 99

```

i

```

Error: Execute cmd: exit=99
Error: Execute cmd: error=</usr/diskless/crco_data/service-lib.sh: line 241:
rpm: command not found
*** Fatal error: ***

```

```
/usr/diskless/crco_data/crco-install: is not a number
>
Error: Execute cmd: exit=99
Error: Execute cmd: error=</tftpboot/service-lib.sh: line 241: rpm: command
not found
*** Fatal error: ***

/tftpboot/install-services: is not a number
```

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```
during command installation script (post)
Errors:
./Bruker/service-lib.sh: line 241: rpm: command not found
is not a number
/home/sermnadmin/Downloads/Topspin-3.2.7/Topspin3.2PL7/xwinstall.d/modscript
s/flexlm.post: terminated with status 99

Output:
```

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```
*** Error ***
```

```
*
Could not show firewall rules for Diskless:
(command: /bin/bash /usr/diskless/prog/bin/setfirewall.sh show)

/usr/diskless/prog/bin/service-lib.sh: line 241: rpm: command not found
Fatal error: is not a number
/usr/diskless/prog/bin/setfirewall.sh: terminated with status 3
```

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```
sermadmin@cie-55-31:~$ topspin
TopSpin 3.5 pl 6 - Copyright (C) 2016 Bruker BioSpin GmbH
Installation directory: "/opt/topspin3.5pl6"
Waiting for FLEXlm license

Cannot checkout TopSpin FLEXlm license
Invalid license file syntax
Feature:      TOPSPIN3
License path: /usr/local/flexlm/Bruker/licenses/license.dat
FLEXlm error: -2,413
For further information, refer to the FLEXlm End User Manual,
available at "www.macrovision.com".
```

The FLEXlm host ID of this machine is 0010182755de

If your workstation controls the spectrometer,
you may run the program without license.
This requires the acquisition server to be ready.

```
checking connection to spectrometer (may last several minutes)
CPR : 2017-05-19 17:30:51.941 +0200
TopSpin is currently running under
the spectrometer service license.
Please apply to receive a correct
license for your workstation
using the online license form:
https://www.bruker.com/nmr_license_requests.html
```

```
error message cprserver: GetSystemPrinterList: dlopen(libcups.so) failed
Exception in thread "main" java.lang.UnsatisfiedLinkError:
/opt/topspin3.5pl6/jre/lib/i386/libawt_xawt.so: libXtst.so.6: cannot open
shared object file: No such file or directory
at java.lang.ClassLoader$NativeLibrary.load(Native Method)
```

```
at java.lang.ClassLoader.loadLibrary0(ClassLoader.java:1941)
at java.lang.ClassLoader.loadLibrary(ClassLoader.java:1824)
at java.lang.Runtime.load0(Runtime.java:809)
at java.lang.System.load(System.java:1086)
at java.lang.ClassLoader$NativeLibrary.load(Native Method)
at java.lang.ClassLoader.loadLibrary0(ClassLoader.java:1941)
at java.lang.ClassLoader.loadLibrary(ClassLoader.java:1845)
at java.lang.Runtime.loadLibrary0(Runtime.java:870)
at java.lang.System.loadLibrary(System.java:1122)
at java.awt.Toolkit$3.run(Toolkit.java:1636)
at java.awt.Toolkit$3.run(Toolkit.java:1634)
at java.security.AccessController.doPrivileged(Native Method)
at java.awt.Toolkit.loadLibraries(Toolkit.java:1633)
at java.awt.Toolkit.<clinit>(Toolkit.java:1668)
at java.awt.Color.<clinit>(Color.java:275)
at de.bruker.nmr.mfw.base.BSplashScreen.<clinit>(BSplashScreen.java:44)
at de.bruker.nmr.mfw.base.MfwBuilder.programMain(MfwBuilder.java:420)
at bruker.bio.start.topspin.Start.main(Start.java:47)
```

premature Java Virtual Machine termination
Program is exiting ...

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Ara el TopSpin ja s'executa,

```

sermnadmin@cie-55-31:~$ topspin
TopSpin 3.5 pl 6 - Copyright (C) 2016 Bruker BioSpin GmbH
Installation directory: "/opt/topspin3.5pl6"
Waiting for FLEXlm license

Cannot checkout TopSpin FLEXlm license
Invalid license file syntax

```

```
Feature: TOPSPIN3
License path: /usr/local/flexlm/Bruker/licenses/license.dat
FLEXlm error: -2,413
For further information, refer to the FLEXlm End User Manual,
available at "www.macrovision.com".
```

The FLEXlm host ID of this machine is 0010182755de

If your workstation controls the spectrometer,
you may run the program without license.
This requires the acquisition server to be ready.

```
checking connection to spectrometer (may last several minutes)
CPR : 2017-05-19 17:33:57.165 +0200
TopSpin is currently running under
the spectrometer service license.
Please apply to receive a correct
license for your workstation
using the online license form:
https://www.bruker.com/nmr_license_requests.html
```

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IPSO: connected to spectrometer subnet

- TCP/IP address = 149.236.99.90
- IPS0 type: AQS
- Tctrl : 1
- Fctrls: 2
- Gctrl1: with digital preemphasis
- Rctrl : 1
- FREDs : none
- DPP : 1

DRU1: AQS DRU Z100977/00582 ECL 04.00

- TCP/IP address = 149.236.99.89
- Firmware Version = 151014
- DRU1 controls AQS-Rack and HPPR/2

AQS: connected to 149.236.99.89:/dev/tty10

Slot	SBSB	Board				ECL	Name	Description
Number	Addr	Type	HW-VS	FW-VS	ID			
2	0x10	0x72	0x3	AW	R	3.0	REC-1	AQS RXAD400 Z102116/705 ECL 03.00
3	0x34	0x2	0x1		X	5.0	REF-1	REF-400 Reference Board for AQS Receiver
4	0x24	0x11	0x8	AE	S	2.0	SGU-1	AQS SGU/2 400 Z103080/00351 ECL 02.00 FPGA-Vs=AF
5	0x25	0x11	0x8	AE	S	2.0	SGU-2	AQS SGU/2 400 Z103080/00340 ECL 02.00 FPGA-Vs=AF
9	0x36	0x5	0		A	6.0	BLA-1	BLA2BB 150/60 20-400 W1345049/1817 ECL 6
21	0x21	0xcf	0		P	1.0	PSD-1	PSD/3 Power Supply Distribution Board
--	0x20	0xd1	0		B	2.0	MASTER	AQS/3 Chassis
1	0x20	0x7	0		B		MASTER	AQS Rack Master

Logical Receiver/SGU configuration:

Router: 1 AQS-Minirouter

Amplifiers (AQS configuration):

RO	Amplifier	Module	Nucleus	Power/W	Switchbox	Name	PN/SN/FW
1	1	1	B	150	N	BLA2BB	150/60 20-400 W1345049/1817/-
2	1	2	B	60	N	BLA2BB	150/60 20-400 W1345049/1817/-
3	2	1	2H	20	N	2HTX in BSMS	-/-/-

Gradient amplifiers at the spectrometer subnet:

BGA1: BGA_W1213891_0040

- TCP/IP address = 149.236.99.93
- Firmware version = 20080206
- Web version = 1.1
- Current limits = 0.0/X, 0.0/Y, -10.0/Z (in A)

BSMS: BSMS/2 connected to ethernet

- TCP/IP address = 149.236.99.20
- ELCB firmware version = 20150423
- SCB channels = 20
- Shim System = BOSS1-SB
- 2HTX Amplifier: available
- Lock: installed

MAS2 control unit: connected to 149.236.99.90:/dev/tty05

- Boot Software : 20031216
- Application Software : 20120329
- Insert Time : 10 s
- Eject Time : 10 s
- Spin Lock Tolerance : 5 Hz

Eurotherm VTU: device connected to 149.236.99.90:/dev/tty04

Preamplifiers :

HPPR2: - HPPR/2 preamplifier connected to 149.236.99.89:/dev/tty10

Type : HPPR/2

Controller: Cover/1 with Application firmware = AQ
15 LEDs for tuning, 15 LEDs for matching

Module 1 : 1H/2H 1H (virtual 50 Ohm reference: 274.9%/85.9deg)
PN=Z003959, SN=00418 from 20070319

Module 2 : 1H/2H 2H (virtual 50 Ohm reference: 116.6%/-71.1deg)
PN=Z003959, SN=00418 from 20070319

Module 3 : XBB19F 2HS
PN=Z003504, SN=00415 from 20070425

Gradient Controller cable connections

GCtrl1: connected via DPP to GREAT 1-10 Gradient Amplifier Unit
W1213891/0040

RF cable connections (assumed)

SGU1 NORM output -> input 2 of transmitter 1 (BLA2BB 150/60 20-400
W1345049/1817)

SGU1 AUX output -> input 1 of transmitter 2 (2HTX in BSMS)

SGU2 NORM output -> input 3 of transmitter 1 (BLA2BB 150/60 20-400
W1345049/1817)

SGU2 AUX output -> TUNE signal input of HPPR

Blanking cable connections (assumed)

transmitter 1 = BLA2BB 150/60 20-400 W1345049/1817:

- amplifier B-150W uses blanking 1
- amplifier B-60W uses blanking 2

```
transmitter 2 = 2HTX in BSMS:
- amplifier 2H-20W uses blanking 3
```

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```
atma_ProbeCorba_impl::start exception: Error:can't stat the wobble process
[NucBB-NucBB- ]
```

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-  `<TOPSPIN_HOME>/exp/stan/nmr/par/users/`
-  `<TOPSPIN_HOME>/exp/stan/nmr/par/teo*` (i `pau*`, `mir*`. Tots es copien a `<TOPSPIN_HOME>/exp/stan/nmr/par/users/`)
-  `<TOPSPIN_HOME>/exp/stan/nmr/lists/pp/users/`
-  `<TOPSPIN_HOME>/exp/stan/nmr/lists/gp/users/` (carpeta buida)
-  `<TOPSPIN_HOME>/exp/stan/nmr/lists/vc/` (carpeta buida)
-  `<TOPSPIN_HOME>/exp/stan/nmr/lists/vd/`
-  `<TOPSPIN_HOME>/exp/stan/nmr/lists/vt/` (carpeta buida)

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```
# 20160422 - MCE - Hide nmrsu user from the login screen
IncludeAll = false
Include = sermnadmin,sermnuab
```

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```
# 20160422 - MCE - Prevent nmrsu user from log-in via ssh
DenyUsers nmrsu
```

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```
root@cie-55-31:/home/sermnadmin/Documents/old-old-hdd-recovery/sdb4-old-old-hdd-opt/topspin/data# rsync -Hvax --progress --stats ./*/opt/topspin3.5pl6/data/
[...]
Number of files: 7,573 (reg: 6,405, dir: 1,168)
Number of created files: 7,573 (reg: 6,405, dir: 1,168)
Number of deleted files: 0
Number of regular files transferred: 6,405
Total file size: 459,101,587 bytes
Total transferred file size: 459,101,587 bytes
Literal data: 459,101,587 bytes
Matched data: 0 bytes
File list size: 65,536
File list generation time: 0.001 seconds
File list transfer time: 0.000 seconds
Total bytes sent: 459,608,752
Total bytes received: 128,851

sent 459,608,752 bytes  received 128,851 bytes  54,086,776.82 bytes/sec
total size is 459,101,587  speedup is 1.00
```

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```
root@cie-55-31:/home/sermnadmin/Documents/old-old-hdd-recovery/sdb4-old-old-hdd-opt/sermn-uab# rsync -Hvax --progress --stats ./data /home/sermnuab/
```

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```
root@cie-55-31:/etc# cat /etc/exports
# /etc/exports: the access control list for filesystems which may be
exported
#
#           to NFS clients.  See exports(5).
#
# Example for NFSv2 and NFSv3:
# /srv/homes          hostname1(rw,sync,no_subtree_check)
hostname2(ro,sync,no_subtree_check)
#
# Example for NFSv4:
```

```
# /srv/nfs4          gss/krb5i(rw,sync,fsid=0,crossmnt,no_subtree_check)
# /srv/nfs4/homes   gss/krb5i(rw,sync,no_subtree_check)
#

# SERMN - EXPORTED FOLDERS

# Export data folder to SeRMN data server
/home/sermnuab/data
sermn.uab.es(ro,sync,no_subtree_check,root_squash,all_squash,insecure)

# Export TopSpin 3.5 pl6 lists folder to SeRMN data server
/opt/topspin3.5pl6/exp/stan/nmr/lists
sermn.uab.es(ro,sync,no_subtree_check,root_squash,all_squash,insecure)

# DO NOT MAKE CHANGES BELOW THIS LINE
/usr/diskless/dl_usr          spect(sync,ro)
/usr/diskless/clients/spect   spect(sync,rw,no_root_squash)
```

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```
root@cie-55-31:/etc# exportfs -av
exportfs: /etc/exports [7]: Neither 'subtree_check' or 'no_subtree_check'
specified for export "spect:/usr/diskless/dl
_usr".
  Assuming default behaviour ('no_subtree_check').
  NOTE: this default has changed since nfs-utils version 1.0.x

exportfs: /etc/exports [8]: Neither 'subtree_check' or 'no_subtree_check'
specified for export "spect:/usr/diskless/cl
ients/spect".
  Assuming default behaviour ('no_subtree_check').
  NOTE: this default has changed since nfs-utils version 1.0.x

exporting spect:/usr/diskless/clients/spect
exporting spect:/usr/diskless/dl_usr
exporting sermn.uab.es:/opt/topspin3.5pl6/exp/stan/nmr/lists
exporting sermn.uab.es:/home/sermnuab/data
```

Llavors provo a muntar les carpetes al servidor de dades i no em dona cap error, i es pot accedir al contingut de les carpetes.

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Nuclei	Probe: 1H-BB/19F BBO+ATMA			Console: SGU & Amp
	ATMA / ATMM (Tune & Match)	Power Limit	PROSOL (Pulses)	CORTAB (RF Calibration)
1H	✓	✓	✓	✓
13C	✓	✓	✓	✓
31P	✓	✓	✓	✓
19F	✓	✓	✓	✓
2H	✓	✓	✓	✓

From: <https://sermn.uab.cat/wiki/> - **SeRMN Wiki**

Permanent link: https://sermn.uab.cat/wiki/doku.php?id=informatica:hp_xw4400_400sb_topspin35&rev=1496408623

Last update: **2017/06/02 15:03**