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The recommended way to install the required TopSpin runtime environment for CentOS-5.x86_64 is:

```
rpm --import ftp://ftp.bruker.de/pub/nmr/CentOS/RPM-GPG-KEY-Bruker-Addon-support
wget ftp://ftp.bruker.de/pub/nmr/CentOS/5/Bruker-Addon/bruker-addon-latest.x86_64.rpm
yum localinstall bruker-addon-latest.x86_64.rpm
```

This will automatically install the configuration and GPG keys that are required for safely installing Bruker AddOn packages.

As default, the repo is disabled. You can either enable it globally, or use `--enablerepo=bruker` when calling yum or the shortcut yumb. Afterwards you may install Bruker specific packages.

For **Processing Workstation**:

```
yum install --enablerepo=bruker bruker-topspin-environ
```

and for **Aquisition**:

```
yum install --enablerepo=bruker bruker-acquisition-environ
```

without the version stamp and you will get always the latest version.

This is a meta-package which installs the TopSpin dependent libraries via its require list:

A full list you get with

```
rpm --requires -qp
'ftp://ftp.bruker.de/pub/nmr/CentOS/5/Bruker-Addon/bruker-topspin-environ-1.0-11.el5.brk.x86_64.rpm'
```

The bruker-acquisition-environ dependencies:

```
rpm --requires -qp ftp://ftp.bruker.de/pub/nmr/CentOS/5/Bruker-  
Addon/bruker-acquisition-environ-1.0-2.el5.brk.noarch.rpm
```

```
/bin/sh  
bootparamd  
bruker-lockdown-session  
bruker-topspin-environ  
dhcp  
nfs-utils  
rpmllib(CompressedFileNames) <= 3.0.4-1  
rpmllib(PayloadFilesHavePrefix) <= 4.0-1  
tftp  
tftp-server  
xinetd
```

All mentioned libraries may add more dependencies during installation

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- Install Bruker software including TopSpin 3.2 and CMC-assist 2.0
- Install GPL licensed GNU GCC 4.5.3 compiler (recommended)

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Com a *setup type* trio

- Data processing, acquisition and automation

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```
/bin/sh /opt/Bruker/mics/setfirewall.sh install
```

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```
Errors from /tftpboot/install-services for IpsoOS:
Installing and configuring DHCP and TFTP service ...
Note: Using 'rpm' for package investigation.
Note: Using init directory /etc/rc.d/init.d
Package dhcp is installed.
Enabling service dhcpd ... [OK]
/etc/dhcpd.conf is not an original Bruker file => saved as
/etc/dhcpd.conf_save ...
/etc/sysconfig/dhcpd is not an original Bruker file => saved as
/etc/sysconfig/dhcpd_save ...
Starting dhcpd: [FAILED]
Package tftp-server is installed.
Enabling service tftp ... [OK]
Installing and configuring DHCP and TFTP service ... [FAILED]
```

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```
[root@cie-55-32 log]# more /etc/sysconfig/networking/devices/ifcfg-
```

```
eth1_spect
ONBOOT=yes
BOOTPROTO=none
NETMASK=255.255.255.0
IPADDR=149.236.99.1
DEVICE=eth1
TYPE=Ethernet
USERCTL=no
IPV6INIT=no
PEERDNS=yes
```

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```
Flexlm will use the TCP/IP prot 1753
If you require a different port number for any reason,
you will have to change the port number in the files
```

```
/usr/local/flexlm/Bruker/licenses/license.dat
and
/etc/services
```

After changing the port number the license manager must be restarted.

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 **incloure captura de pantalla**

```
*** Notice ***

Error: Execute cmd: exit=1

Error: Execure cmd: error=<installing and configuring DHCP and TFTP service
...
Note: Using 'rpm' for package investigation.
Note: Using init directory /etc/rc.d/init.d
Package dhcp is installed
Enabling service dhcpd ... [OK]
/etc/dhcpd.conf is not an original Bruker file => saved as
/etc/init.d/dhcpd.conf_save ...
/etc/sysconfig/dhcpd is not an original Bruker file => saved as
/etc/sysconfig/dhcpd_save ...
```

```
Starting dhcpd: [FAILED]
Package tftp-server is installed
Enabling service tftp ... [OK]
Installing and configuring DHCP and TFTP service ... [FAILED]
```

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Please close existing windows or logout and login again before you launch the program. </note>

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```
>> u

Autoboot: Waiting to load
```

```
bfs()/usr/diskless/clients/spect/root/unix.r4600_std (CTRL-C to abort)
loading
```

```
Fast ethernet (FEN): 100 MBit/s, half duplex
Obtaining /usr/diskless/clients/spect/root/unix.r4600_std from server
cie-55-32
799504+98976+804288 entry: 0x80021000
CPU: MIPS R4700 Processor Chip Revision: 1.0
FPU: MIPS R4700 Floating Point Unit [CP1] Revision: 1.0
```

```
RISC/os Release 4_52 mips Version UMIPS
Total real memory  = 16777216
Available memory   = 14442496
Root on nfs file :,Swap on nfs file :
Fast ethernet (FEN) : 100 MBit/s, full duplex
whoami RPC call failed with status 3
PANIC: vfs_mountroot: cannot mount root
syncing disks... done
```

```
Fast ethernet (FEN): 100 MBit/s, half duplex
```

```
BRUKER CCU R4600 Monitor
Version 5.40 MIPS OPT Thu May 17 12:18:17 MSZ 2001 mar
Memory size: 16777216 (0x1000000) bytes, 16 MB
Icache size: 16384 (0x4000) bytes
Dcache size: 16384 (0x4000) bytes
>>
```

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```
whoami RPC call failed with status 3

PANIC: vfs_mountroot: cannot mount root
syncing disks... done
```

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```
[root@cie-55-32 etc]# more /etc/bootparams
spect    root=ASP_ST2:/usr/diskless/clients/spect/root
swap=ASP_ST2:/usr/diskless/clients/spect/swapfile
dump=ASP_ST2:/usr/diskless/clients/spect/dump
[root@cie-55-32 etc]# ps ax | grep boot
 6306 pts/2    R+      0:00 grep boot
[root@cie-55-32 etc]# /etc/init.d/bootparamd
Usage: /etc/init.d/bootparamd {start|stop|status|restart|reload}
[root@cie-55-32 etc]# /etc/init.d/bootparamd status
rpc.bootparamd is stopped
```

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```
rpc status 15 RPC: Program not registered
```

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```
[root@cie-55-32 etc]# /etc/init.d/nfs status
rpc.mountd is stopped
nfsd is stopped
rpc.rquotad is stopped
```

i per resoldre el problema torno a engegar el servidor NFS,

```
[root@cie-55-32 etc]# /etc/init.d/nfs start
Starting NFS services:                                [ OK ]
Starting NFS quotas:                                [ OK ]
Starting NFS daemon:                                [ OK ]
Starting NFS mountd:                                [ OK ]
Stopping RPC idmapd:                                [ OK ]
Starting RPC idmapd:                                [ OK ]
```

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```
Root on nfs file :,Swap on nfs file :  
Fast ethernet (FEN) : 100 MBit/s, full duplex  
hostname: spect  
domainname: (none)  
Root fstype nfs  
Available memory    = 12795904  
  
The system is coming up.  Please wait.  
ASP_ST2:/usr/diskless/clients/spect/var mounted on /var  
ASP_ST2:/usr/diskless/dl_usr mounted on /usr  
Fast ethernet (FEN) : 100 MBit/s, full duplex  
Internet daemons: portmap inetd.  
NFS daemons: biod lockd statd.  
The system is ready.  
  
spect Console login:
```

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DRU1 Firmware requires upgrade

✔ fet!

```
Firmware information for DRU1:  
  Current version on DRU1 : 2008-05-09  
  Minimum required version: 2013-01-25  
  Available for update    : 2013-01-25  
  
Please proceed as follows:  
1) Wait until 'cf' is finished!  
2) Start Mozilla Web browser and enter the URL 'http://149.236.99.89'  
3) Click 'AQS Firmware Setup' -> 'Load new DRU firmware'  
4) Browse to '/opt/topspin3.2/conf/instr/servtool/dru'  
5) Select the file 'dru_firmware_130125.gz and install the firmware.  
6) Before starting any acquisition you must execute 'ii' in TopSpin!  
  
=====
```

11 August 2015 17:02:55.174 +0200
Topspin Version = 3.2 (of May 14 2013),build 1791
JVM Version = 1.7.0_13 Oracle Corporation
JVM Total memory = 91 MB

JVM Free memory = 20 MB

AQS_RX1 Firmware requires upgrade

✔ fet!

ATTENTION: the firmware on AQS_RX1 is too old and must be updated!

FPGA Firmware information for Receiver1:

Current version on Receiver 1 : AT

Minimum required version: AW

Available for update : AW

Please proceed as follows:

- 1) Open a terminal shell and start 'UniTool'
- 2) Hit the <Enter> key or enter 'aq's' as device name
- 3) Enter 16 as SBSB address
- 4) Read the warning, then hit the <Enter> key
- 5) Select 'Receiver Service Menu' and hit the <Enter> key
- 6) Select 'Manual Download' and hit the <Enter> key
- 7) Enter 'birds/rxs_aw.hex' and hit the <Enter> key
- Wait until the update is finished - this may take up to 30 min.
- 8) Enter 'X' repeatedly until UniTool terminates

If the update is not done 'wobb' may fail with RGP_RX error.

```
=====
11 August 2015 17:03:06.600 +0200
Topspin Version = 3.2 (of May 14 2013),build 1791
JVM Version      = 1.7.0_13 Oracle Corporation
JVM Total memory = 91 MB
JVM Free  memory = 19 MB
```

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```
Application Version : AW
Application Date    : 12-10-29
Application Checksum : 0x038
Board Type         : 0x72
Board Hardware Version : 2
Board ECL          : 1.1
Boot Version       : AD
Boot Date          : 03-08-28
```

MAS Unit configuration

✔ fet!

```
MAS control unit: Firmware 981102
Your MAS2 is configured for compatibility mode.
To enable support of all features, please reconfigure
to native mode and do 'cf' again.
```

```
=====
11 August 2015 17:08:27.000 +0200
Topspin Version = 3.2 (of May 14 2013),build 1791
JVM Version      = 1.7.0_13 Oracle Corporation
JVM Total memory = 93 MB
JVM Free  memory = 51 MB
```

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```
Cannot open device for Preemphasis Unit:
Opening preemphasis unit failed:
Initial open of tty (/dev/tty06) without delay failed: No such device or
address
Probably the GREAT is switched off
```

```
=====
11 August 2015 17:08:32.050 +0200
Topspin Version = 3.2 (of May 14 2013),build 1791
JVM Version      = 1.7.0_13 Oracle Corporation
JVM Total memory = 93 MB
JVM Free  memory = 50 MB
```

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Warning:

```
Current MAS firmware "20051025" does not support MAS2!  
You must update to firmware "20080601" or newer to use the features of the  
MAS2  
For a firmware update open the MAS-display, switch to the config panel,  
enable the edit mode and press the download button.
```

```
=====  
12 August 2015 10:25:46.032 +0200  
Topspin Version = 3.2 (of May 14 2013), build 1791  
JVM Version      = 1.7.0_13 Oracle Corporation  
JVM Total memory = 90 MB  
JVM Free  memory = 41 MB
```

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```
*****
*                >>> bsms service toolbox <<<                *
:
                >>> Main Menu <<<
0  Keyboard On/Off
1  Init BSMS
2  Check/Download all boards
3  Show all versions
:
F  Service...
Q Quit      ? Help      ! System      Your choice: 2

Checking out need for download...
:
D O W N L O A D   -   R E P O R T

Keyboard is up to date
SCB-M is up to date
SCB-R is up to date
Lock needs autodownload
SLCB needs autodownload
CPU is up to date
communicationtest (15 seconds) ...

Estimated autodownload time: 0 h 33 min
Press 'y' to start autodownload. Any other key aborts. y

Lock ( 40/ 40/ 16) autodownload with firmware (lockap.hex)
Reading file 'lockap.hex'.
check HEX-file (Intel) ...
Erasing FLASH memory...

Programming record: 4969 of 4970
Download done.
54% done
SLCB ( 72/ 40/ -1) autodownload with firmware (slcbad.hex)
Reading file 'slcbad.hex'.
check HEX-file (Intel) ...
Erasing FLASH memory...

Programming record: 4123 of 4123
Download done.
100% done

BSMS Init

2. download-check
Boards in system: CPLSTUD
:
D O W N L O A D   -   R E P O R T
```

```
Keyboard is up to date
SCB-M is up to date
SCB-R is up to date
Lock is up to date
SLCB is up to date
CPU is up to date
```

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```
Eurotherm VTU: device connected to spect:/dev/tty05
```

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```
Connecting SBSB address 80 (0x50).
HPPR/2 Main Menu
=====
[0] Init HPPR/2
[1] Show Version
:
[X] Exit
    your choice: 1
HPPR/2 Cover Module Version Data
```

```
-----  
Type      Date      V0      V1      CheckSum  
App SW    071018      A0      -      00EA  
Boot SW    040928      AD      -      00A6  
Hardware   -          1      0      -
```

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HPPR/2 Main Menu

```
=====
[0] Init HPPR/2
[1] Show Version
[A] Auto Download
```

:

```
[X] Exit
your choice: A
```

Check for latest firmware and instructions:

```
ftp://ftp.bruker.ch/pub/NMR/download/servtools/firmware/hppr/
Press '.' to abort, <ENTER> to start download:
```

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Logical Receiver/SGU configuration:

Router: 1 Avance-Router

Amplifiers (default configuration):

R0	Amplifier	Module	Nucleus	Power/W	Switchbox	Name	PN/SN/FW
1	1	1	X	100	N	(I=00) BLAXH	50 -/-/-
2	1	2	1H	50	N	(-""-) BLAXH	50 -/-/-

- Type: BSMS/1
- absolute lockshift frequency = changed from 11426996 Hz to 11428531 Hz Hz
- GAB type = GAB
- GAB current limits = 0.0/X, 0.0/Y, 10.0/Z (in A)

Lock: device connected to spect:/dev/tty03

MAS2 control unit: connected to spect:/dev/tty07

- Boot Software : 20031216
- Application Software : 20130613
- Insert Time : 10 s
- Eject Time : 10 s
- Spin Lock Tolerance : 5 Hz

Eurotherm VTU: device connected to spect:/dev/tty05

Preamplifiers :

HPPR: - HPPR/2 preamplifier connected to spect:/dev/tty10

Type : HPPR/2

Controller: Application firmware = AQ

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MAS error

The MAS Firmware does not allow to define the maximum values for DemandedRate, DrivePressure, and BearingPressure, therefore the 'edhead' parameters

PHMROT=15000

PHMDRV=0

PHMBRG=0

are ignored

Please define the 'installed Probe' in the 'masdisp' window to set the correct maximum values and spin up procedure for the installed MAS Probe!

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```
at de.bruker.cpr.CPRClientPOA._invoke(CPRClientPOA.java:43)
at org.openorb.adapter.poa.POA.dispatch(Unknown Source)
at org.openorb.net.AbstractServerRequest.dispatch(Unknown Source)
at org.openorb.net.ServerManagerImpl.serve_request(Unknown Source)
at org.openorb.net.ServerManagerImpl.thread_pool_main(Unknown Source)
at org.openorb.net.ServerManagerImpl.access$200(Unknown Source)
at org.openorb.net.ServerManagerImpl$PoolThread.run(Unknown Source)
```

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```
In file included from /opt/topspin3.2/gnu/lib/gcc/i686-pc-linux-
gnu/4.5.3/include-fixed/features.h:328:0,
                 from /usr/include/stdio.h:28,
                 from nmr_save.c:30:
/usr/include/gnu/stubs.h:7:27: fatal error: gnu/stubs-32.h: No such file or
directory
compilation terminated.
exit 1
```

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```
In file included from /opt/topspin3.2/gnu/lib/gcc/i686-pc-linux-
gnu/4.5.3/include-fixed/features.h:328:0,
                 from /usr/include/stdio.h:28,
                 from 2dflshift.c:30:
/usr/include/gnu/stubs.h:7:27: fatal error: gnu/stubs-32.h: No such file or
directory
compilation terminated.
exit 1
```

La cerca a Google de “stubs-32.h centos 5” em proporciona la solució,

Fix for the RHEL/CentOS 5.x for GCC gnu/stub-32.h missing error Type the following yum command:

```
# yum -y install glibc-devel.i386
```

[x86_64 Linux Error: gnu/stub-32.h Missing Error and Solution](#)

La provo

```
[root@cie-55-32 gnu]# yum -y install glibc-devel.i386
```

Loaded plugins: fastestmirror, security

Loading mirror speeds from cached hostfile

- * base: mirror.uv.es
- * epel: mirror.uv.es
- * extras: mirror.uv.es
- * rpmforge: fr2.rpmfind.net
- * rpmforge-extras: fr2.rpmfind.net
- * updates: mirror.uv.es

Setting up Install Process

Resolving Dependencies

--> Running transaction check

---> Package glibc-devel.i386 0:2.5-123.el5_11.1 set to be updated

--> Finished Dependency Resolution

Dependencies Resolved

=====			
=====			
Package	Arch	Version	
Repository	Size		
=====			
=====			
Installing:			
glibc-devel	i386	2.5-123.el5_11.1	updates
2.1 M			

Transaction Summary

=====

Install	1 Package(s)
Upgrade	0 Package(s)

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```
protocol version mismatch - is your shell clean?  
(see the rsync man page for an explanation)  
rsync error: protocol incompatibility (code 2) at compat.c(62)  
rsync: connection unexpectedly closed (0 bytes received so far) [Receiver]  
rsync error: error in rsync protocol data stream (code 12) at io.c(226)  
[Receiver=3.1.1]
```

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```
sent 1,740,752 bytes  received 6,170,330,889 bytes  6,095,873.23 bytes/sec  
total size is 6,162,532,109  speedup is 1.00
```

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```
sent 118,960 bytes  received 242,854,959 bytes  8,835,415.24 bytes/sec
total size is 242,305,330  speedup is 1.00
```

```
$ rsync -aHxv --protocol=26 sermn@192.168.2.7:/home/sermn .
sermn@192.168.2.7's password:
receiving file list ... done
sermn/
sermn/.adobe/
sermn/.adobe/Acrobat/
:
sermn/Desktop/USB Stick
sermn/Desktop/www.redhat.com
sermn/teo/esquire6000.pdf
```

```
sent 4,832 bytes  received 15,339,082 bytes  2,789,802.55 bytes/sec
total size is 15,314,312  speedup is 1.00
```

```
$ rsync -aHxv --protocol=26 teo@192.168.2.7:/home/teo .
teo@192.168.2.7's password:
receiving file list ... done
teo/
teo/.adobe/
teo/.adobe/Acrobat/
teo/.adobe/Acrobat/7.0/
:
teo/NMRSIM_SESSION/wave/rect.shp
teo/NMRSIM_SESSION/wave/sinc5neg.shp
```

```
sent 4,288 bytes  received 2,123,317 bytes  607,887.14 bytes/sec
total size is 2,101,812  speedup is 0.99
```

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where RO is the router output number.

- *TSHOME/ntp/<nucleus>/FCU<no>/<ampl. descr.>/<date>/*
 - *Amp_Power_res* - theoretical and experimental amplifier output voltage
 - *Amp_Phase_res* - experimental amplifier phase values
 - *Att_verif pha.<expno>* - amplifier verification result
 - *Elec_Hard_Testfile* - all Cortab results of the current date

.conf* - hardware configurations (output of Hardware Configuration) *functions. - Cortab functions (output of Hardware Configuration)

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```
90 degree pulse import for
/opt/topspin3.1/conf/instr/spect/prosol/33/generic/19F_dec failed!
```

```

/opt/topspin3.1/conf/instr/spect/prosol/33/generic/1H_obs succeeded
/opt/topspin3.1/conf/instr/spect/prosol/33/generic/1H_dec succeeded
/opt/topspin3.1/conf/instr/spect/prosol/33/generic/2H@o_obs succeeded
/opt/topspin3.1/conf/instr/spect/prosol/33/generic/2H@o_dec succeeded
/opt/topspin3.1/conf/instr/spect/prosol/33/generic/11B_obs succeeded
/opt/topspin3.1/conf/instr/spect/prosol/33/generic/11B_dec succeeded
/opt/topspin3.1/conf/instr/spect/prosol/33/generic/13C_obs succeeded
/opt/topspin3.1/conf/instr/spect/prosol/33/generic/13C_dec succeeded
/opt/topspin3.1/conf/instr/spect/prosol/33/generic/15N_obs succeeded
/opt/topspin3.1/conf/instr/spect/prosol/33/generic/19F_obs succeeded
/opt/topspin3.1/conf/instr/spect/prosol/33/generic/29Si_obs succeeded
/opt/topspin3.1/conf/instr/spect/prosol/33/generic/29Si_dec succeeded
/opt/topspin3.1/conf/instr/spect/prosol/33/generic/31P_obs succeeded
/opt/topspin3.1/conf/instr/spect/prosol/33/generic/31P_dec succeeded

```

Get more information about failed imports by typing "hist"!

Please check the succeeded imports experimentally!

The customized import from within edpsol provides more chances.

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- <TOPSPIN_HOME>/exp/stan/nmr/par/users/
- <TOPSPIN_HOME>/exp/stan/nmr/lists/pp/users/
- <TOPSPIN_HOME>/exp/stan/nmr/lists/gp/users/
- <TOPSPIN_HOME>/exp/stan/nmr/lists/vc/
- <TOPSPIN_HOME>/exp/stan/nmr/lists/vd/
- <TOPSPIN_HOME>/exp/stan/nmr/lists/vt/
- <TOPSPIN_HOME>/exp/stan/nmr/lists/wave/users/

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PXELINUX 2.11 2004-08-16 Copyright (C) 1994-2004 H. Peter Anvin