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- Select the type of acquisition: High Resolution Systems
- Available spectrometer or datastation configurations
 - Select the configuration you want to use: spect
- Select the items you want to install
 - Install Pulse Programs
 - Install Bruker AU Programs
 - [unchecked] Recompile All User AU Programs
 - Install Library CPD Programs
 - Install Library Gradient Files
 - Install Library Shape Files
 - Convert Standard Parameter Sets
 - Install Standard Scaling Region Files
 - Install Bruker Python Programs
- Select the basic frequency of your spectrometer
 - Basic Frequency: 400.13
- Select the pre-scan delay DE
 - Default pre-scan-delay (us): 6.5
- Select the plotter paper format
 - Paper format: A4 / Letter

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- SAVE INSTALLATION FILES
 - Location of backup file: /opt/topspin3.6.3/nmr_backup * Installation to be saved (TopSpin home) : /opt/topspin3.6.3
 - Spectrometer configuration (e.g. spect): spect
 - Display default information: checked
 - Display additional information: unchecked
- Execute “Save installation specific files” periodically: [Automatic Backup]

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```
retrieving path to diskless (may take some time)...
found /usr/diskless
prepare backup: Copy /usr/diskless/prog/logfiles to
/opt/topspin3.6.3/nmr_backup/tmp_backup_dir/logfiles
extract and backup /opt/topspin3.6.3/conf/instr/servtool/bsmstool/boss/
=>
extract and backup /opt/topspin3.6.3/conf/instr/servtool/bsmstool/boss/wb/
now extracting from dir
/opt/topspin3.6.3/conf/instr/servtool/bsmstool/boss/wb/...
[...]...
Every file of the following directories will be saved:

/opt/topspin3.6.3/QTP/
/opt/topspin3.6.3/conf/global/
/opt/topspin3.6.3/conf/instr/...
[...]
/opt/topspin3.6.3/db/user/
/opt/topspin3.6.3/data/final/nmr/protocolfiles/

Additionally the following files will be saved (if existing):

/opt/topspin3.6.3/conf/instr/spect/uxnmr.par
/opt/topspin3.6.3/conf/instr/spect/uxnmr.info
[...]
/opt/topspin3.6.3/prog/server/export.conf

-----
start to generate back-up copy...

Now writing zip file
/opt/topspin3.6.3/nmr_backup/completeDirectories.zip1...
[...]
Backing up of NMR data was successful!
The backup file

/opt/topspin3.6.3/nmr_backup/nmr_backup_20210824-1726.zip

has been generated!

#####
```

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- SAVE USER FILES
 - Location of backup file: /opt/topspin3.6.3/nmr_backup
 - Overwrite existing backup file unchecked
 - Display default information: checked
 - Display additional information: unchecked
- Execute Save user files periodically:

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```
retrieving path to diskless (may take some time)...
found /usr/diskless
```

```
Every non-Bruker file in /home/sermnadmin/.topspin-cie-56-18.uab.cat/ will
be saved:
```

```
Every non-Bruker file in /opt/topspin3.6.3/exp/stan/nmr/au/src/ will be
saved:
```

```
now extracting from dir /opt/topspin3.6.3/exp/stan/nmr/au/src/...
```

```
Now writing zip file
/opt/topspin3.6.3/nmr_backup/extractedDirectories.zip1...
```

```
Every non-Bruker file in /opt/topspin3.6.3/exp/stan/nmr/lists/base_info/
will be saved:
```

```
now extracting from dir /opt/topspin3.6.3/exp/stan/nmr/lists/base_info/...
```

```
[...]
```

```
Every non-Bruker file in /opt/topspin3.6.3/exp/stan/nmr/lists/vd/ will be
saved:
```

```
now extracting from dir /opt/topspin3.6.3/exp/stan/nmr/lists/vd/...
```

```
Every non-Bruker file in /opt/topspin3.6.3/exp/stan/nmr/lists/vp/ will be
saved:
```

```
now extracting from dir /opt/topspin3.6.3/exp/stan/nmr/lists/vp/...
```

Every non-Bruker file in /opt/topspin3.6.3/exp/stan/nmr/lists/vt/ will be saved:

```
now extracting from dir /opt/topspin3.6.3/exp/stan/nmr/lists/vt/...
```

Every non-Bruker file in /opt/topspin3.6.3/exp/stan/nmr/lists/wave/ will be saved:

```
now extracting from dir /opt/topspin3.6.3/exp/stan/nmr/lists/wave/...
```

Now writing zip file

```
/opt/topspin3.6.3/nmr_backup/extractedDirectories.zip7...
```

Every non-Bruker file in /opt/topspin3.6.3/exp/stan/nmr/par/ will be saved:

```
now extracting from dir /opt/topspin3.6.3/exp/stan/nmr/par/...
```

Now writing zip file

```
/opt/topspin3.6.3/nmr_backup/extractedDirectories.zip8...
```

Every non-Bruker file in /opt/topspin3.6.3/exp/stan/nmr/py/ will be saved:

```
now extracting from dir /opt/topspin3.6.3/exp/stan/nmr/py/...
```

Now writing zip file

```
/opt/topspin3.6.3/nmr_backup/extractedDirectories.zip9...
```

Now writing zip file

```
/opt/topspin3.6.3/nmr_backup/nmr_backup_sermnadmin_20210812-1137.zip...
```

```
#####
```

Backing up of NMR data was successful!

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exp/stan/nmr/lists/ (multiples carpetes)

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