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```
$ sh ./topspin-3.6.4-linux.sh
Verifying archive integrity... All good.
Uncompressing TopSpin 3.6.4 100%
NOTE: using DISPLAY=:0
Please enter root password if prompted for it.
Password:
NOTE: using DISPLAY=:0
NOTE: using TCL_LIBRARY=/tmp/selfgz1602632396/linux/tcl-8.6.9/lib/tcl8.6
NOTE: using TK_LIBRARY=/tmp/selfgz1602632396/linux/tk-8.6.9/lib/tk8.6
NOTE: running /tmp/selfgz1602632396/linux/tk-8.6.9/bin/wish8.6 -f
/tmp/selfgz1602632396/xwinstall.d/swim/lib/xwinstall.tcl --
NOTE: using log file: /tmp/install.log
```

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- Install Bruker NMR Software
- Install GPL licensed GNU GCC 4.9.1 compiler (recommended)

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- Data processing, acquisition and automation  
Installs Data Processing, Acquisition and Automation Components for Spectroscopy

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- Select a spectrometer configuration to import
  - spect in /opt/topspin3.6.3

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```
*** Warning ***

Cannot import the selected spectrometer configuration (got invalid
```

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```
The IPSO computer controlling the spectrometer has to be rebooted for the installed Diskless package to take effect.
```

```
Reboot IPSO computer automatically? [Yes]
```

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```
Error: *** Error ***
```

```
Installation of Spectrometer Operating System: Could not reboot the IPSO
computer:
connect failed connection refused
```

You have to reset or reboot the IPSO computer at IP address 149.236.99.90 manually before the installed Diskless package takes effect.

To do this, either

- press the red button on the IPSO computer, or
- start TopSpin after this installation has been completed, type 'ha', open the IPSO web site, navigate to 'Administration' and click on Reboot.

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```
# systemctl status tftp
● tftp.service - Tftp Server
   Loaded: loaded (/usr/lib/systemd/system/tftp.service; indirect; vendor
  preset: disabled)
   Active: inactive (dead)
     Docs: man:in.tftpd

# systemctl enable tftp
Created symlink from /etc/systemd/system/sockets.target.wants/tftp.socket to
/usr/lib/systemd/system/tftp.socket.

# systemctl start tftp

# systemctl status tftp
● tftp.service - Tftp Server
   Loaded: loaded (/usr/lib/systemd/system/tftp.service; indirect; vendor
  preset: disabled)
   Active: active (running) since Thu 2021-08-26 16:02:06 CEST; 2s ago
     Docs: man:in.tftpd
  Main PID: 13576 (in.tftpd)
   CGroup: /system.slice/tftp.service
           └─13576 /usr/sbin/in.tftpd -s /var/lib/tftpboot
```

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```
2021-08-26 15:57:11.616 +0200
Cannot checkout TopSpin CodeMeter license
No matching license entry found for feature "TOPSPIN3" mapped to product
code 1011013000 .
CmContainer Entry not found, Error 200.
```

```
CodeMeter Runtime Version: 7.10.4196
empty CodeMeter Server Search List Configuration
```

```
If your workstation controls the spectrometer,
you may run the program without license.
This requires the acquisition server to be ready.
```

```
checking connection to spectrometer (may last several minutes)
CPR : 2021-08-26 15:57:13.104 +0200
```

TopSpin is currently running under the spectrometer service license. Please apply to receive a correct license for your workstation using the online license form: [https://www.bruker.com/nmr\\_license\\_requests.html](https://www.bruker.com/nmr_license_requests.html)

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- Spectrometer Description
  - Description: Avance III 400 @SeRMN
- Spectrometer Data
  - 1H Spectrometer frequency: 400.130 MHz
- Security Options
  - Enable power check: yes (option checked)
- Optional Standard Devices
  - MAS2 control unit: connected to: tty05
  - Eurotherm VTU: device connected to: tty04

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TopSpin 3.6.4 - Copyright (C) 2021 Bruker BioSpin GmbH