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```
$ sh ./topspin-3.6.4-linux.sh
Verifying archive integrity... All good.
Uncompressing TopSpin 3.6.4 100%
NOTE: using DISPLAY=:0
Please enter root password if prompted for it.
Password:
NOTE: using DISPLAY=:0
NOTE: using TCL_LIBRARY=/tmp/selfgz1602632396/linux/tcl-8.6.9/lib/tcl8.6
NOTE: using TK_LIBRARY=/tmp/selfgz1602632396/linux/tk-8.6.9/lib/tk8.6
NOTE: running /tmp/selfgz1602632396/linux/tk-8.6.9/bin/wish8.6 -f
/tmp/selfgz1602632396/xwinstall.d/swim/lib/xwinstall.tcl --
NOTE: using log file: /tmp/install.log
```

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- Install Bruker NMR Software
- Install GPL licensed GNU GCC 4.9.1 compiler (recommended)

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- Data processing, acquisition and automation  
Installs Data Processing, Acquisition and Automation Components for Spectroscopy

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- Select a spectrometer configuration to import
  - spect in /opt/topspin3.6.3

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```
*** Warning ***

Cannot import the selected spectrometer configuration (got invalid
```

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```
The IPSO computer controlling the spectrometer has to be rebooted for the installed Diskless package to take effect.
```

```
Reboot IPSO computer automatically? [Yes]
```

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```
Error: *** Error ***
```

```
Installation of Spectrometer Operating System: Could not reboot the IPSO
computer:
connect failed connection refused
```

You have to reset or reboot the IPSO computer at IP address 149.236.99.90 manually before the installed Diskless package takes effect.

To do this, either

- press the red button on the IPSO computer, or
- start TopSpin after this installation has been completed, type 'ha', open the IPSO web site, navigate to 'Administration' and click on Reboot.

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```
# systemctl status tftp
● tftp.service - Tftp Server
   Loaded: loaded (/usr/lib/systemd/system/tftp.service; indirect; vendor
  preset: disabled)
   Active: inactive (dead)
     Docs: man:in.tftpd

# systemctl enable tftp
Created symlink from /etc/systemd/system/sockets.target.wants/tftp.socket to
/usr/lib/systemd/system/tftp.socket.

# systemctl start tftp

# systemctl status tftp
● tftp.service - Tftp Server
   Loaded: loaded (/usr/lib/systemd/system/tftp.service; indirect; vendor
  preset: disabled)
   Active: active (running) since Thu 2021-08-26 16:02:06 CEST; 2s ago
     Docs: man:in.tftpd
  Main PID: 13576 (in.tftpd)
   CGroup: /system.slice/tftp.service
           └─13576 /usr/sbin/in.tftpd -s /var/lib/tftpboot
```

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```
++ test -z ''
++ export 'JAVA_ARGS=-
splash:/opt/topspin3.6.4/classes/prop/loadinggraphic_600.png -Xms64M -
Xmx200M -XX:+HeapDumpOnOutOfMemoryError -Xverify:none -
DsecondVirtualMachine=n -DXWINMRHOME="/opt/topspin3.6.4" -
DCOMPUTERNAME=cie-54-214 -DUSERHOME_DOT_TOPSPIN="/home/sermnadmin/.topspin-
cie-54-214" -Duser.language= -Dawt.useSystemAAFontSettings=lcd -
Djava.endorsed.dirs="/opt/topspin3.6.4/classes/lib/jacorb" -
DTOPSPIN_COMM_DEBUG=0 -Djava.util.Arrays.useLegacyMergeSort=true -jar
"/opt/topspin3.6.4/classes/lib/topspin.jar"
++ JAVA_ARGS='-splash:/opt/topspin3.6.4/classes/prop/loadinggraphic_600.png
-Xms64M -Xmx200M -XX:+HeapDumpOnOutOfMemoryError -Xverify:none -
DsecondVirtualMachine=n -DXWINMRHOME="/opt/topspin3.6.4" -
DCOMPUTERNAME=cie-54-214 -DUSERHOME_DOT_TOPSPIN="/home/sermnadmin/.topspin-
cie-54-214" -Duser.language= -Dawt.useSystemAAFontSettings=lcd -
Djava.endorsed.dirs="/opt/topspin3.6.4/classes/lib/jacorb" -
DTOPSPIN_COMM_DEBUG=0 -Djava.util.Arrays.useLegacyMergeSort=true -jar
"/opt/topspin3.6.4/classes/lib/topspin.jar"
+ '[' '' -z '' ']'
+ StartProgram /opt/topspin3.6.4/prog/mod/cpr -c demo
+ '[' /opt/topspin3.6.4/prog/mod/cpr = topspin ']'
+ case "$machine" in
+ export LD_LIBRARY_PATH=/opt/topspin3.6.4/prog/shlib
+ LD_LIBRARY_PATH=/opt/topspin3.6.4/prog/shlib
+ Program=/opt/topspin3.6.4/prog/mod/cpr
+ shift
++ dirname /opt/topspin3.6.4/prog/mod/cpr
+ d=/opt/topspin3.6.4/prog/mod
+ '[' /opt/topspin3.6.4/prog/mod = . ']'
+ '[' '' -f /opt/topspin3.6.4/prog/mod/cpr ']'
+ '[' '' -x /opt/topspin3.6.4/prog/mod/cpr ']'
+ MyExec /opt/topspin3.6.4/prog/mod/cpr -c demo
+ exec /opt/topspin3.6.4/prog/mod/cpr -c demo
TopSpin 3.6.4 - Copyright (C) 2021 Bruker BioSpin GmbH
Installation directory: "/opt/topspin3.6.4"
Waiting for CodeMeter license

2021-08-26 15:57:11.616 +0200
Cannot checkout TopSpin CodeMeter license
No matching license entry found for feature "TOPSPIN3" mapped to product
code 1011013000 .
CmContainer Entry not found, Error 200.

CodeMeter Runtime Version: 7.10.4196
empty CodeMeter Server Search List Configuration

If your workstation controls the spectrometer,
you may run the program without license.
This requires the acquisition server to be ready.

checking connection to spectrometer (may last several minutes)
```

```
CPR : 2021-08-26 15:57:13.104 +0200
TopSpin is currently running under
the spectrometer service license.
Please apply to receive a correct
license for your workstation
using the online license form:
https://www.bruker.com/nmr_license_requests.html
```

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- Spectrometer Description
  - Description: Avance III 400 @SeRMN
- Spectrometer Data
  - 1H Spectrometer frequency: 400.130 MHz
- Security Options
  - Enable power check: yes (option checked)
- Optional Standard Devices
  - MAS2 control unit: connected to: tty05
  - Eurotherm VTU: device connected to: tty04

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