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```
$ sh ./topspin-3.6.4-linux.sh
Verifying archive integrity... All good.
Uncompressing TopSpin 3.6.4  100%
NOTE: using DISPLAY=:0
Please enter root password if prompted for it.
Password:
NOTE: using DISPLAY=:0
NOTE: using TCL_LIBRARY=/tmp/selfgz1602632396/linux/tcl-8.6.9/lib/tcl8.6
NOTE: using TK_LIBRARY=/tmp/selfgz1602632396/linux/tk-8.6.9/lib/tk8.6
NOTE: running /tmp/selfgz1602632396/linux/tk-8.6.9/bin/wish8.6 -f
/tmp/selfgz1602632396/xwinstall.d/swim/lib/xwinstall.tcl --
NOTE: using log file: /tmp/install.log
```

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- Install Bruker NMR Software
- Install GPL licensed GNU GCC 4.9.1 compiler (recommended)

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- Data processing, acquisition and automation  
Installs Data Processing, Acquisition and Automation Components for Spectroscopy

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- Select a spectrometer configuration to import
  - spect in /opt/topspin3.6.3

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\*\*\* Warning \*\*\*

Cannot import the selected spectrometer configuration (got invalid

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```
The IPSO computer controlling the spectrometer has to be rebooted for the installed Diskless package to take effect.
```

```
Reboot IPSO computer automatically? [Yes]
```

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```
Error: *** Error ***
```

Installation of Spectrometer Operating System: Could not reboot the IPSO computer:  
connect failed connection refused

You have to reset or reboot the IPSO computer at IP address 149.236.99.90 manually before the installed Diskless package takes effect.

To do this, either

- press the red button on the IPSO computer, or
- start TopSpin after this installation has been completed, type 'ha', open the IPSO web site, navigate to 'Administration' and click on Reboot.

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```
# systemctl status tftp
● tftp.service - Tftp Server
  Loaded: loaded (/usr/lib/systemd/system/tftp.service; indirect; vendor
  preset: disabled)
    Active: inactive (dead)
      Docs: man:in.tftpd

# systemctl enable tftp
Created symlink from /etc/systemd/system/sockets.target.wants/tftp.socket to
/usr/lib/systemd/system/tftp.socket.

# systemctl start tftp

# systemctl status tftp
● tftp.service - Tftp Server
  Loaded: loaded (/usr/lib/systemd/system/tftp.service; indirect; vendor
  preset: disabled)
    Active: active (running) since Thu 2021-08-26 16:02:06 CEST; 2s ago
      Docs: man:in.tftpd
    Main PID: 13576 (in.tftpd)
   CGroup: /system.slice/tftp.service
           └─13576 /usr/sbin/in.tftpd -s /var/lib/tftpboot
```

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2021-08-26 15:57:11.616 +0200  
Cannot checkout TopSpin CodeMeter license

```
No matching license entry found for feature "TOPSPIN3" mapped to product code 1011013000 .
CmContainer Entry not found, Error 200.
```

```
CodeMeter Runtime Version: 7.10.4196
empty CodeMeter Server Search List Configuration
```

If your workstation controls the spectrometer,  
you may run the program without license.  
This requires the acquisition server to be ready.

```
checking connection to spectrometer (may last several minutes)
```

```
CPR : 2021-08-26 15:57:13.104 +0200
```

TopSpin is currently running under  
the spectrometer service license.  
Please apply to receive a correct  
license for your workstation  
using the online license form:  
[https://www.bruker.com/nmr\\_license\\_requests.html](https://www.bruker.com/nmr_license_requests.html)

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- Spectrometer Description
  - Description: Avance III 400 MHz SB - SeRMN, UAB
- Spectrometer Data
  - 1H Spectrometer frequency: 400.130 MHz
- Security Options
  - Enable power check: yes (option checked)
- Optional Standard Devices
  - MAS2 control unit: connected to: tty05
  - Eurotherm VTU: device connected to: tty04

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 Error setting parameter "LIFTMODE": Parameter does not exist  
(Original message = Error setting parameter "LIFTMODE": Parameter does not exist)

=====

26 August 2021 16:17:50.630 +0200

Topspin Version = TopSpin 3.6.4 (of 2021-07-26 14:56:34), build 95

OS Name = Linux 3.10.0-1160.25.1.el7.x86\_64  
JVM Version = 1.8.0\_242 (64 bit) AdoptOpenJDK  
JVM Total memory = 236 MB  
JVM Free memory = 142 MB

```
java.lang.Exception: Stack trace
    at
bruker.bio.topspin.bfw.gui.warn.AbstractMessage.getDetailsString(A
bstractMessage.java:875)
    at
bruker.bio.topspin.bfw.gui.warn.AbstractMessage.setMsgSource(Absr
actMessage.java:837)
    at
bruker.bio.topspin.bfw.gui.warn.AbstractMessage.<init>(AbstractMes
age.java:821)
    at de.bruker.nmr.mfw.base.BWarning.<init>(BWarning.java:59)
    at de.bruker.nmr.mfw.base.BError.<init>(BError.java:50)
    at de.bruker.nmr.sc.util.Misc.makeDialog(Misc.java:119)
    at de.bruker.nmr.sc.util.Misc.showMsg(Misc.java:98)
    at de.bruker.nmr.sc.util.Misc.showMsg(Misc.java:72)
    at de.bruker.nmr.sc.util.Misc.showMsg(Misc.java:59)
    at
de.bruker.nmr.sc.bsms.impl.BsmsCorbaUtil$1.run(BsmsCorbaUtil.java:
147)
    at java.lang.Thread.run(Thread.java:748)
```

From:  
<https://sermn.uab.cat/wiki/> - SeRMN Wiki

Permanent link:  
[https://sermn.uab.cat/wiki/doku.php?id=informatica:hp\\_z230\\_400sb\\_topspin364\\_1&rev=1629987724](https://sermn.uab.cat/wiki/doku.php?id=informatica:hp_z230_400sb_topspin364_1&rev=1629987724)

Last update: 2021/08/26 16:22

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TopSpin 3.6.4 - Copyright (C) 2021 Bruker BioSpin GmbH