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```
$ sh ./topspin-3.6.3-linux.sh
Verifying archive integrity... All good.
Uncompressing TopSpin 3.6.3 100%
NOTE: using DISPLAY=:0.0
Please enter root password if prompted for it.
Password: *****
NOTE: using DISPLAY=:0.0
NOTE: using TCL_LIBRARY=/tmp/selfgz158608332/linux/tcl-8.6.9/lib/tcl8.6
NOTE: using TK_LIBRARY=/tmp/selfgz158608332/linux/tk-8.6.9/lib/tk8.6
NOTE: running /tmp/selfgz158608332/linux/tk-8.6.9/bin/wish8.6 -f
/tmp/selfgz158608332/xwinstall.d/swim/lib/xwinstall.tcl --
NOTE: using log file: /tmp/install.log
```

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- Install Bruker NMR Software
- Install GPL licensed GNU GCC 4.9.1 compiler (recommended)

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- Data processing, acquisition and automation
Installs Data Processing, Acquisition and Automation Components for Spectroscopy

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- Spectrometer Description
 - Description: Avance III 600 @SeRMN
- Spectrometer Data
 - 1H Spectrometer frequency: 600.130 MHz
- Security Options
 - Enable power check: option not checked
- Optional Standard Devices
 - Bruker Automatic Changer (BACS) connected to: ttty01
 - Variable Temperature Unit connected to: 149.236.99.20
- Options for the Sample Changer
 - Sample lift controlled by: Sample Changer
 - Delay between insert of sample and next command: 10 s
 - Start experiment while sample changer in motion: option not checked

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BSMS Error:

SampleTransporter: Emergency button pressed.

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