



HyperSense®

DNP for NMR



HyperSense®
Delivering Sensitivity to NMR

The Business of Science™



HyperSense[®]

Delivering Sensitivity to NMR

HyperSense - the *in-vitro* DNP Polariser


HyperSense allows the user to drive NMR sensitivity to levels never seen before. The sensitivity improvements attainable with ¹³C, ¹⁵N and ²⁹Si have the potential to increase the breadth of your Nuclear Magnetic Resonance (NMR) experiments and your applications.

Examples are:

- Shorter data collection times and/or the use of much smaller amounts of sample material.
- Benefits for quaternary carbons that will improve structure confirmation and elucidation. This could provide more reliable NMR spectral assignments in a shorter time.
- The ability using direct detection of ¹³C to identify complex biomixtures in metabonomics.
- Using the enhanced signal to carry out chemical kinetic experiments.

HyperSense, an instrument that configures easily with your laboratory and existing spectrometer.


¹³C-DNP of acetaminophen natural abundance, 2 μmol (300 μg)



What is DNP?


Dynamic Nuclear Polarisation (DNP) is the method for obtaining strongly polarised nuclear spins in solution, thereby delivering sensitivity enhancement for your application.

In a glassy state the free electron of the radical, at low temperature, in the presence of a magnetic field can transfer the polarisation to nearby nuclei using microwave irradiation.



Enhanced ¹³C Spectral Content

- After polarisation, carbon nuclei are enhanced. This makes detection quicker compared to conventional NMR.



4 hours polarisation, single scan
Limit of detection S/N of 10 on the smallest peak F

Working with **HyperSense**[®] in your Laboratory

1. Preparing your sample for DNP

Dissolve your sample in a glassing agent (e.g. DMSO/H₂O). Add proprietary radical (typically 15mM). Place sample in polariser sample cup (10 - 200µl). Load sample into polariser.



2. Sample Conditioning

Sample in polarising region. Sample cooled to $\geq 1.4K$

3. Polarising

Irradiate unpaired electrons with microwaves (over 15 minutes to 4 hours) to transfer polarisation

4. Transfer



Sample dissolution using hot solvent. Polarised sample is transferred

5. Acquisition

HyperSense triggers NMR acquisition



The Automated Process of Polarisation





What does **HyperSense**[®] deliver?

Examples of challenging direct detection

¹⁵N Direct Detection

¹⁵N of mixture of L-alanine t-BOC, glycine t-BOC and N-acetyl glycine (10mg of each)



- Single scan
- 5 hours polarisation



²⁹Si Direct Detection

²⁹Si of Volasil-244 (9.5mg)

- Single scan
- 5 hours polarisation



Features of HyperSense:

High levels of sample polarisation

- Enables the user to directly detect ¹³C, ¹⁵N and ²⁹Si more quickly than with conventional NMR
- Requires smaller amounts of sample, down to µg level

After polarisation some carbons are particularly enhanced

- Enhanced spectral content improving peak assignment

Separate hyperpolarisation module


- Enables dual functionality of conventional NMR and DNP-assisted spectroscopy

Case Studies using HyperSense®



1. Chemical Kinetics

Low sensitivity limits kinetic studies using NMR. Long acquisition times restrict most studies to ^1H . For the first time HyperSense provides real time multiple ^{13}C spectra from a single sample. Direct ^{13}C detection improves spectral resolution and makes the observation of carbon reaction centres possible.

Relative amount of product




From ratio of peak integral, at 172.3ppm (reaction:control) kinetic rates can be determined¹



- $^{13}\text{C}_1$ Ethylpyruvate polarised for 100min
- Sample transferred using methanol to NMR magnet
- 5mm NMR tube contains 75mM NaOH solution
- ^{13}C spectra acquired one per second

2. Detection of biologically relevant molecules in mixtures



NMR of complex mixtures can be difficult due to low sensitivity and spectral crowding for ^1H spectra. Conventional ^{13}C NMR offers greater chemical shift dispersion but with lower sensitivity.


The single-scan ^{13}C spectrum of a compound mixture (citrate, acetamide, tartrate, vitamins B1 and B3) shows the work being carried out at Birmingham University into DNP NMR of complex mixtures². Using HyperSense the NMR spectrum of multiple molecules can be enhanced, so aiding interpretation.

Case Studies of HyperSense®

3. Enhanced 2D NMR

2D NMR provides information rich correlated spectra. The single enhanced spectrum produced by HyperSense can be used in conjunction with the SNR requirements of the ultrafast 2D single scan method developed by Professor Lucio Frydman^{3,4}.

The figure shows a ^1H - ^{15}N HSQC spectrum of ^{15}N labelled urea, 130 μM . Natural abundance



materials at low concentration can also be used to record both ^1H - ^{13}C and ^1H - ^{15}N HSQC spectra⁵.


4. Research into DNP

DNP is a new and exciting area of NMR spectroscopy which in itself provides huge scope for research.

Sample polarisation is affected by: temperature; microwave power; radical identity and concentration; polarisation time and the nature of the sample.

By way of example the figure below shows the effect of temperature. With increasing temperature there is an inverse dependence on polarisation and the rate of polarisation build-up.

HyperSense offers the user the ability to change the polarisation parameters⁶.



Other applications include:

- Structure verification
- Impurity identification
- Ligand screening
- Biomarker discovery

Summary

HyperSense will transform your research using ^{13}C , ^{15}N and ^{29}Si . It will enable you to surpass previous sensitivity barriers and increase the applicability of NMR spectroscopy by:

- Increasing signal to noise ratio up to 10,000 times⁷

- Reducing experimental runtimes for ^{13}C and ^{15}N applications

- Dramatically lowering limit of signal detection

- Switching rapidly and simply between NMR and DNP NMR



Let us work with you to show you the benefits of DNP NMR with HyperSense. At Oxford Instruments we have several demonstration facilities across the world staffed with expert scientists to enable you to assess the benefits of DNP NMR to your applications.

Performance Specification and Services

System performance

Nuclei available:	^{13}C , ^{15}N , ^{29}Si and other spin $1/2$ nuclei
Sample polarisation times:	Typically 15 minutes to 6 hours
NMR spectrometer compatibility:	All current spectrometer platforms
NMR RF probe requirements:	Operation with 5 and 10mm probes
Dissolution solvents:	Water and methanol
Polarising agent:	'Trityl radical' typical usage 0.2-5mg per sample

Instrument specification

Microwave source:	Integrated 94GHz with 0.5GHz sweep, (user selectable) up to 100mW
Magnetic field strength:	3.35T actively shielded
Sample temperature:	Selectable $<3.9\text{K}$, $\geq 1.4\text{K}^{\circ}$
Helium refill volume:	65L from minimum level
Nitrogen refill volume:	60L from minimum level
Helium consumption during dissolutions:	Typical $<2\text{L}$ per dissolution ⁹
Minimum operating ceiling height:	$<3.0\text{m}$
System weight including cryogenics:	600kg approx
Trigger output for NMR spectrometer:	5V signal

Services

Compressed air/nitrogen:	6 bar min
Helium gas:	6 bar min, grade 99.999% pure from cylinder
Exhaust line:	Vacuum exhaust line or suitable safety compliant filter
Electrical:	Europe: 230V, 50Hz, 2.5A single phase and 415V, 50Hz 20A three phase USA: 110V, 60Hz, 3A single phase and 208V, 60Hz, 20A three phase Japan: 100V, 50-60Hz, 3A single phase and 200V, 60Hz, 20A three phase

Dimensions

Footprint:	W 1668mm (~66") x D 1038mm (~41")
Height:	2860mm (~113")



Oxford Instruments Molecular Biotoools

UK

Tubney Woods, Abingdon
Oxon OX13 5QX
Tel: +44 (0) 1865 393 200
Fax: +44 (0) 1865 393 333
E-mail: molecularbiotoools@oxinst.co.uk

Japan

2-11-6, Tomioka, Koto-ku
Tokyo 135-0047
Tel: +81 (0) 3 5245 3261
Fax: +81 (0) 3 5245 4472
E-mail: molecularbiotoools@oxinst.co.jp

USA

8403 Cross Park Drive, Suite 3F
Austin Texas 78754
Tel: +1 512 339 0640
Fax: +1 512 339 0620
E-mail: molecularbiotoools@ma.oxinst.com

click onto www.oxford-instruments.com for more information

References and notes

- (1) Detection of Minor Isomers and Study of Chemical Kinetics by DNP NMR – Presented at SMASH 2006, Damir Blazina, Steven Reynolds, Oxford Instruments Molecular Biotoools Ltd., Oxon, UK.*
 - (2) Probing DNP Enhancement for Biological Samples - Presented at ENC 2006, Christian Ludwig¹, Abdul Hamid-Emwas¹, Damir Blazina², Steven Reynolds², Andrew Sowerby¹, Ulrich Gunter¹. ¹CR UK Institute for Cancer Studies, University of Birmingham, Birmingham, UK, ²Oxford Instruments Molecular Biotoools Ltd., Oxon, UK.*
 - (3) Frydman, L.; Scherf, T.; Lupulescu, A. *Proc. Nat. Acad. Sci.* **2002**, *99*, 15858.
 - (4) Using the 2D single scan technology requires a license from the Weizmann Institute of Science.
 - (5) *Ultrafast Two-Dimensional Nuclear Magnetic Resonance Spectroscopy of Hyperpolarised Solutions*, L. Frydman and D. Blazina, *Nature Physics*, in press.
 - (6) Characterising Solid-State DNP, Rob Slade, Graham Hutton, Damir Blazina.*
 - (7) Ardenkjaer-Larsen, J.H.; Fridlung, B.; Gram, A.; Hansson, G.; Hansson, L.; Lerche, M. H.; Servin, R.; Thaning, M.; Golman, K. *Proc. Nat. Acad. Sci.* **2003**, *100*, 10158.
 - (8) As measured by instrument control system with 65% helium insert level at microwave power 100mW over 30 minutes.
 - (9) Averaged over 10 repeatable runs with a temperature of 1.4K over a 6 hour polarisation and dissolution, 100mW microwave power with helium level at 65%.
- *Please see www.oxford-instruments.com.

Ref: Hyp/0307

This publication is the copyright of Oxford Instruments Molecular Biotoools Ltd and provides outline information only, which (unless agreed by the company in writing) may not be used, applied or reproduced for any purpose or form part of any order or contract or regarded as the representation relating to the products or services concerned. Oxford Instruments' policy is one of continued improvement. The company reserves the right to alter, without notice the specification, design or conditions of supply of any product or service. HyperSense is a registered trademark of Oxford Instruments Molecular Biotoools Ltd. and the Oxford Instruments logo is a trademark of Oxford Instruments plc or its subsidiaries.

HyperSense is not marketed for *in-vivo* use.

www.oxford-instruments.com

