

# Changelog dels canvis de la 1.3.0 a la 1.3.1

Aquests són els canvis fets de la versió 1.3.0 a la versió 1.3.1. Al llistat s'esmenta,

- el nom del fitxer modificat
- la versió del fitxer al [CVS de Sourceforge](#) (**extended:**)
- una descripció breu dels canvis fets

## ChangeLog file

Aquestes són les línies al fitxer *ChangeLog* que descriuen de forma general els canvis entre les 1.3.0 i 1.3.1. A l'apartat següent es poden consultar les modificacions a nivell de fitxer.

- bugfix: **user cannot insert a booking id (in the URL) if the id doesn't exist in database**

## Root files (/)

### index.php

- **Extended: 1.64.1**
- added another *ParseConfig* to update before auth object

## General files (inc/)

### inc/jsfunctions.php

- **Extended: 1.9.1**
- deleted js function 'disp\_confirm()' for language action (because this function will be in menu.php)

### inc/menu.php

- **Extended: 1.35.3**
- added new header (iconNMR) for iconNMR actions
- added js function disp\_confirm() in \_constructLanguageForm()

### inc/permissions.php

- **Extended: 1.3.3**
- added new permissions for iconNMR: BBPERM\_ADMIN\_QUEUES (to administrate queues) and BBPERM\_SAMPLES (to manage samples of the queues)

## BB files (inc/bb/)

### inc/bb/auth.php

- **Extended: 1.41.3**
- new function getQueues for iconNMR module
- new functions: isQueueAdmin() (like isInstrumentAdmin()) and manageSamplesPermitted() for iconNMR module

### inc/bb/configreader.php

- **Extended: 1.12.2**
- BUMBLEBEEVERSION variable is taken from bumblebee.ini
- added \$this→iconNMR to verify if iconNMR module is on

### inc/bb/conservation.php

- Initial revision
- Conservation object for iconNMR

### inc/bb/dangerousness.php

- Initial revision
- Dangerousness object for iconNMR

### inc/bb/experiment.php

- Initial revision
- Experiment object for iconNMR

### inc/bb/instrument.php

- **Extended: 1.32**
- list of users with access to instrument doesn't display if we are creating an instrument
- fix 'emailonbooking' checkbox

### inc/bb/queue.php

- Initial revision
- Queue object for iconNMR module
- When you create a queue you also have to select the experiments and solvents of it

## **inc/bb/sendsample.php**

- Initial revision of send sample object (iconNMR module)

## **inc/bb/settings.php**

- **Extended: 1.7.1**
- added 'Extra modules' section (with iconNMR module)

## **inc/bb/solvent.php**

- Initial revision
- Solvent object for iconNMR module

## **inc/bb/user.php**

- **Extended: 1.32.3**
- added queue form options for iconNMR module (similar to instrument form)

## **Actions files (inc/actions/)**

### **inc/actions/actions.php**

- **Extended: 1.25.3**
- added conditionals for iconNMR module
- added actions: sendsample, my samples, manage samples, download samples file, edit queues, edit solvents, edit experiments, edit conservations and edit dangerousness (for iconNMR module)

### **inc/actions/actionfactory.php**

- **Extended: 1.30.2**
- added require\_once for actions of iconNMR module

### **inc/actions/calendar.php**

- **Extended: 1.1.3**
- if instrid is invalid it shows the error using reportAction function

### **inc/actions/conservations.php**

- Initial revision
- Conservations action for iconNMR module

## **inc/actions/dangerousness.php**

- Initial revision
- Dangerousness action for iconNMR module

## **inc/actions/experiments.php**

- Initial revision
- Experiment action for iconNMR module

## **inc/actions/mysamples.php**

- Initial revision
- To view a list of your sended/processed samples like “My bookings”
- If you are an admin queue, you will be able to view deleted samples too

## **inc/actions/password.php**

- **Extended: 1.15.1**
- added <h2> title

## **inc/actions/queues.php**

- Initial revision
- Queue action for iconNMR module

## **inc/actions/samples.php**

- Initial revision of samples form (iconNMR module)
- to manage samples and create iconNMR file

## **inc/actions/samples.php**

- Initial revision of samplesfile form (iconNMR module)
- to download iconNMR file

## **inc/actions/sendsample.php**

- Initial revision of send sample form (iconNMR module)
- with a powerful function: `_MultiChoiceJoin()` (could be an object in the future?)

## inc/actions/solvents.php

- Initial revision
- Solvent action for iconNMR module

## inc/actions/view.php

- **Extended: 1.62.3**
- added selectQueue() function for iconNMR module

## Formlib files (inc/formslib/)

### inc/formslib/checkboxtablelist.php

- **Extended: 1.12.1**
- add variable \$linkAll: boolean make all columns in the table link to the target URL

### inc/formslib/nondbrow.php

- **Extended: 1.13**
- style changes in displayInTable function

## Configuration files (config/)

### config/bumblebee.ini

- **Extended: 1.65.3**
- added 'version' parameter (for BUMBLEBEEVERSION)
- added new section [modules]
- added new module: iconNMR

## Install files (install/)

### install/setup-tables.sql

- **Extended: 1.11.3**
- new tables: iconNMR\_queues and iconNMR\_permissions
- new tables: iconNMR\_solvents, iconNMR\_experiments, iconNMR\_dangerousness and iconNMR\_conservations
- new tables: iconNMR\_queue\_experiment and iconNMR\_queue\_solvent to make a link between queues-solvents and queues-experiments
- new tables: iconNMR\_samples (sample form)

- new tables: iconNMR\_sample\_experiment, iconNMR\_sample\_dangerousness and iconNMR\_sample\_conservation to make a link between them.

From:  
<https://sermn.uab.cat/wiki/> - **SeRMN Wiki**

Permanent link:  
[https://sermn.uab.cat/wiki/doku.php?id=sermn\\_wiki:userpages:marta:changelog\\_130\\_a\\_131&rev=1337612026](https://sermn.uab.cat/wiki/doku.php?id=sermn_wiki:userpages:marta:changelog_130_a_131&rev=1337612026)

Last update: **2012/05/21 16:53**

