

Changelog dels canvis de la 1.3.0 a la 1.3.1

Aquests són els canvis fets de la versió 1.3.0 a la versió 1.3.1. Al llistat s'esmenta,

- el nom del fitxer modificat
- la versió del fitxer al [CVS de Sourceforge](#) (**extended:**)
- una descripció breu dels canvis fets

ChangeLog file

Aquestes són les línies al fitxer *ChangeLog* que descriuen de forma general els canvis entre les 1.3.0 i 1.3.1. A l'apartat següent es poden consultar les modificacions a nivell de fitxer.

- feature: iconNMR from Bruker BioSpin module (turn on/off through Configure System option)
- bugfix: user cannot insert a booking id (in the URL) if it doesn't exist in database
- feature: BUMBLEBEEVERSION is taken from bumblebee.ini.
- bugfix: list of users with access to instrument doesn't display if we are creating an instrument
- bugfix: in calendar action, if instrid is invalid will show an error
- feature: some styles change on tables and checkboxtablelist action (link option)
- feature: new section [modules] in bumblebee.ini to create modules like iconNMR.

Root files (/)

index.php

- **Extended: 1.64.1**
- added another *ParseConfig* to update before auth object

General files (inc/)

inc/jsfunctions.php

- **Extended: 1.9.1**
- deleted js function 'disp_confirm()' for language action (because this function will be in menu.php)

inc/menu.php

- **Extended: 1.35.3**
- added new header (iconNMR) for iconNMR actions
- added js function disp_confirm() in _constructLanguageForm()

inc/permissions.php

- **Extended: 1.3.3**
- added new permissions for iconNMR: BBPERM_ADMIN_QUEUES (to administrate queues) and BBPERM_SAMPLES (to manage samples of the queues)

BB files (inc/bb/)

inc/bb/auth.php

- **Extended: 1.41.3**
- new function getQueues for iconNMR module
- new functions: isQueueAdmin() (like isInstrumentAdmin()) and manageSamplesPermitted() for iconNMR module

inc/bb/configreader.php

- **Extended: 1.12.2**
- BUMBLEBEEVERSION variable is taken from bumblebee.ini
- added \$this→iconNMR to verify if iconNMR module is on

inc/bb/conservation.php

- Initial revision
- Conservation object for iconNMR

inc/bb/dangerousness.php

- Initial revision
- Dangerousness object for iconNMR

inc/bb/experiment.php

- Initial revision
- Experiment object for iconNMR

inc/bb/instrument.php

- **Extended: 1.32**
- list of users with access to instrument doesn't display if we are creating an instrument
- fix 'emailonbooking' checkbox

inc/bb/queue.php

- Initial revision
- Queue object for iconNMR module
- When you create a queue you also have to select the experiments and solvents of it

inc/bb/sendsample.php

- Initial revision of send sample object (iconNMR module)

inc/bb/settings.php

- **Extended: 1.7.1**
- added 'Extra modules' section (with iconNMR module)

inc/bb/solvent.php

- Initial revision
- Solvent object for iconNMR module

inc/bb/user.php

- **Extended: 1.32.3**
- added queue form options for iconNMR module (similar to instrument form)

Actions files (inc/actions/)

inc/actions/actions.php

- **Extended: 1.25.3**
- added conditionals for iconNMR module
- added actions: sendsample, my samples, manage samples, download samples file, edit queues, edit solvents, edit experiments, edit conservations and edit dangerousness (for iconNMR module)

inc/actions/actionfactory.php

- **Extended: 1.30.2**
- added require_once for actions of iconNMR module

inc/actions/calendar.php

- **Extended: 1.1.3**

- if instrid is invalid it shows the error using reportAction function

inc/actions/conservations.php

- Initial revision
- Conservations action for iconNMR module

inc/actions/dangerousness.php

- Initial revision
- Dangerousness action for iconNMR module

inc/actions/experiments.php

- Initial revision
- Experiment action for iconNMR module

inc/actions/mysamples.php

- Initial revision
- To view a list of your sended/processed samples like “My bookings”
- If you are an admin queue, you will be able to view deleted samples too

inc/actions/password.php

- **Extended: 1.15.1**
- added <h2> title

inc/actions/queues.php

- Initial revision
- Queue action for iconNMR module

inc/actions/samples.php

- Initial revision of samples form (iconNMR module)
- to manage samples and create iconNMR file

inc/actions/samplesfile.php

- Initial revision of samplesfile form (iconNMR module)
- to download iconNMR file

inc/actions/sendsample.php

- Initial revision of send sample form (iconNMR module)
- with a powerful function: `_MultiChoiceJoin()` (could be an object in the future?)

inc/actions/solvents.php

- Initial revision
- Solvent action for iconNMR module

inc/actions/view.php

- **Extended: 1.62.3**
- added `selectQueue()` function for iconNMR module

Formlib files (inc/formslib/)

inc/formslib/checkboxtablelist.php

- **Extended: 1.12.1**
- add variable `$linkAll`: boolean make all columns in the table link to the target URL

inc/formslib/nondbrow.php

- **Extended: 1.13**
- style changes in `displayInTable` function

Configuration files (config/)

config/bumblebee.ini

- **Extended: 1.65.3**
- added 'version' parameter (for BUMBLEBEEVERSION)
- added new section [modules]
- added new module: iconNMR

Install files (install/)

install/setup-tables.sql

- **Extended: 1.11.3**
- new tables: `iconNMR_queues` and `iconNMR_permissions`

- new tables: iconNMR_solvents, iconNMR_experiments, iconNMR_dangerousness and iconNMR_conservations
- new tables: iconNMR_queue_experiment and iconNMR_queue_solvent to make a link between queues-solvents and queues-experiments
- new tables: iconNMR_samples (sample form), iconNMR_file
- new tables: iconNMR_sample_experiment, iconNMR_sample_dangerousness and iconNMR_sample_conservation to make a link between them.

From:
<https://sermn.uab.cat/wiki/> - **SeRMN Wiki**

Permanent link:
https://sermn.uab.cat/wiki/doku.php?id=sermn_wiki:userpages:marta:changelog_130_a_131&rev=1337616973

Last update: **2012/05/21 18:16**

