Pure Shift NMR Covariance

MERCK

André Fredia, Pau Nolisa, Carlos Cobasb, Gary E. Martinc and Teodor Parella

Universitat Autònoma de Barcelona


^a Servei de Ressonància Magnètica Nuclear, Universitat Autònoma de Barcelona, Barcelona, Catalonia. ^b Mestrelab Research, Santiago de Compostela, Spain. ^c NMR Structure Elucidation. Merck & Co. Inc., Rahway, USA.

Introduction

The development of novel experimental strategies to significantly enhance signal resolution by broadband homodecoupling is a current topic of high interest in ¹H NMR spectroscopy¹. A number of different building blocks have been implemented into 1D and 2D homo- and heteronuclear experiments in order to provide resolution-enhanced pure chemical shift ¹H NMR spectra, where signals appear collapsed to singlets. On the other hand, Covariance processing methods have been used to generate challenging NMR spectral representations². We present here the first attempts towards a general solution to generate *Pure Shift NMR spectra by using Generalized Indirect Covariance* (psGIC) co-processing^{3,4}. The current strategy is based on the calculation of a new 2D psGIC spectrum from the combination of a parent homo- or heteronuclear spectrum and a reference 2D F1-homodecoupled ¹H-¹H correlation spectrum only showing diagonal cross-peaks (DIAG), which share a common ¹H frequency dimension. Using psGIC, the F1 dimension in the DIAG spectrum is transferred to the F2 dimension of the parent spectrum, thus generating a new pure shift 2D spectrum.

Methodology

The key point of this proposal is a highly resolved 2D DIAG spectrum that can be obtained by two methods: i) experimentally, from a modification of the recent 2D PSYCHE-TOCSY pulse sequence with omission of the DIPSI-2 pulse train⁵ (Fig. 1A); ii) a faster method reconstructs it from an experimental 1D broadband homodecoupled ¹H spectrum⁶ followed by the make2D* reconstruction (Mnova) (Fig. 1B). As proof of concept, the conventional 2D HSQC and 2D DIAG spectra of strychnine were combined following the GIC formalism (λ =1). This co-processing directly provides a reconstructed pure shift HSQC (psGIC-HSQC) spectrum showing well defined singlet signals for all cross-peaks (Fig. 2).



Step 2: Covariance **Processing**





Figure 2


Pure-Shift NMR Spectra

This GIC procedure can be used to reconstruct any pure shift spectrum, even for some experiments that can be very challenging or impossible to acquire through experimental pulse schemes. Some examples are illustrated by the psGIC-HSQC-TOCSY (Fig. 3) and psGIC-HSQMBC (Fig. 4) spectra.



Aliphatic region corresponding to (A) the regular 2D HSQC-TOCSY and B) the psGIC-HSQC-TOCSY spectra of strychnine generated by GIC. A comparison of the internal F2 projections at the top of each 2D spectrum and the selected C14 rows in C) shows the absence of J multiplicity and the excellent reproducibility of signal intensities for all cross-peaks.

Figure 3.








Figure 4. Aliphatic region corresponding to (A) the regular HSQMBC and (B) the psGIC-HSQMBC spectra of strychnine generated by GIC. A comparison of the internal F2 projections at the top of each 2D spectrum and the selected C14 rows in (C) shows the absence of coupling pattern and the excellent reproducibility signal intensities for all crosspeaks.

Multiplicity-edited NMR Spectra

The same approach can be applied using the ¹H-¹H 2D pure shift ME-DIAG (ME-psDIAG) spectrum (Fig. 5B) synthetically generated from internal F2-projected (Fig. 5A) spectrum 2D ME-psHSQC experimental spectrum recorded with BIRDhomodecoupling during based acquisition. Thus, ME-GIC versions of the classical COSY (Fig. 5D), TOCSY, HSQC-TOCSY or HMBC spectra can be obtained from GIC combination with the ME-psDIAG 8 7 6 5 4 3 2 1 spectrum.





On the other hand, a ME-psHMBC-COSY could spectrum reconstructed by GIC between the standard HMBC and the transposed ME-psCOSY spectra (Fig. 6), giving a major number of correlations that can be useful to observe missing two-bond correlations in the original HMBC spectrum or to detect up to four- and five bonds away longrange correlations.

- ✓ Simple and fast processing tool.
- ✓ Complement for conventional dataset analysis
- ✓ General application to a broad range of NMR experiments.
- ✓ Generation of novel synthetic pure-shift spectra, even for NMR experiments that are "impossible" nowadays.
- ✓ Access to unfeasible multiplicity-edited COSY, TOCSY or HMBC spectra

References:

¹L. Castañar, T. Parella, Magn. Reson. Chem., 53 (2015) 399. ² M. Jaeger, R.L.E.G. Asper, Ann. Rep. NMR Spectroscopy, 83 (2014) 272.

³ A. Fredi *et al., J. Magn. Reson.*, 266 (2016) 16.

- ⁴ A. Fredi *et al., J. Magn. Reson.*, 270 (2016) 161 ⁵ M. Foroozandeh, et al., J. Am. Chem. Soc., 136 (**2014**) 11867.
 - ⁶ M. Foroozandeh, et al., Angew. Chem. Intl. Ed.., 53 (**2014**) 6990.

Acknowledgements

Financial support for this research provided by MINECO of Spain (project CTQ2015-64436-P) is gratefully acknowledged. We also thank to the Servei de Ressonància Magnètica Nuclear, Universitat Autònoma de Barcelona, for allocating instrument time to this project and to CNPq-Brasil for the PhD scholarship.

Figure 1 make2D is Mnova script (available in Mnova 11) that synthesizes 2D-DIAG spectra by automatically deconvolving a ps-1D spectrum and generating 2D diagonal peaks using Lorentzian-type lineshapes