

# Tranferència de la documentació del BioSpec i PharmaScan

Transfereixo el contingut de la carpeta */opt/Technical Documentation for BioSpec and PharmaScan Systems* de l'ordinador actual a la nova estació de treball,

```
root@cie-50-205:/opt# rsync -Hvax --verbose --progress --stats --
protocol="30" \
-e "ssh" --rsync-path="sudo rsync" sermn@cie-55-30.uab.cat:/opt/Technical* .
| \
tee 20160909-1530-rsync-pv40-docs.log
opening connection using: ssh -l sermn cie-55-30.uab.cat "sudo rsync" --
server --sender -vvlHogDtpxe.iLsfx . "/opt/Technical*" (10 args)
sermn@cie-55-30.uab.cat's password:
receiving file list ...
expand file_list to 4000 bytes, did move
220 files to consider
delta-transmission enabled
Technical Documentation for BioSpec and PharmaScan Systems/
Technical Documentation for BioSpec and PharmaScan Systems/technical_ DOCU/
Technical Documentation for BioSpec and PharmaScan Systems/technical_
DOCU/Amplifier Cabinet (Bruker-Siemens B-GS 2000-500)/
Technical Documentation for BioSpec and PharmaScan Systems/technical_
DOCU/Amplifier Cabinet (Bruker-Siemens B-GS 2000-500)/Active Decoupling/
Technical Documentation for BioSpec and PharmaScan Systems/technical_
DOCU/Amplifier Cabinet (Bruker-Siemens B-GS 2000-500)/Cabinet Wiring/
:
:
Technical Documentation for BioSpec and PharmaScan Systems/technical_
DOCU/RF-Data_Cabinet/Shim_power_supply/T00654_T7012.pdf
      2,016,014 100%    5.72MB/s    0:00:00 (xfr#142, to-chk=2/220)
Technical Documentation for BioSpec and PharmaScan Systems/technical_
DOCU/System_Interconnections/T01733_T9513.pdf
      6,855,183 100%    7.10MB/s    0:00:00 (xfr#143, to-chk=0/220)
total: matches=0 tag_hits=0 false_alarms=0 data=349780731

rsync[17820] (receiver) heap statistics:
arena:          270336    (bytes from sbrk)
ordblks:        4        (chunks not in use)
smlblks:        0
hblks:          3        (chunks from mmap)
hblkhd:        798720    (bytes from mmap)
allmem:        1069056    (bytes from sbrk + mmap)
usmlblks:       0
fsmblks:       0
uordblks:      219840    (bytes used)
fordblks:       50496    (bytes free)
keepcost:      41952    (bytes in releasable chunk)
```

```
rsync[17773] (generator) heap statistics:
arena:          270336   (bytes from sbrk)
ordblks:        3       (chunks not in use)
smlblks:        2
hblks:          2       (chunks from mmap)
hblkhd:         532480   (bytes from mmap)
allmem:         802816   (bytes from sbrk + mmap)
usmlblks:       0
fsmblks:        64
uordblks:       154272   (bytes used)
fordblks:       116064   (bytes free)
keepcost:       107520   (bytes in releasable chunk)
```

Number of files: 220 (reg: 143, dir: 77)  
Number of regular files transferred: 143  
Total file size: 349,780,731 bytes  
Total transferred file size: 349,780,731 bytes  
Literal data: 349,780,731 bytes  
Matched data: 0 bytes  
File list size: 7,281  
Total bytes sent: 2,304  
Total bytes received: 349,836,292

sent 2,304 bytes received 349,836,292 bytes 10,140,249.16 bytes/sec  
total size is 349,780,731 speedup is 1.00

## Verificació

Tipus de fitxer	Ordinador actual	Nou ordinador
Fitxer normal	143	143 ✓
Carpeta	77	77 ✓
Enllaç simbòlic	0	0 ✓

From: <https://sermn.uab.cat/wiki/> - **SeRMN Wiki**

Permanent link: [https://sermn.uab.cat/wiki/doku.php?id=informatica:hp\\_z240mt\\_biospec\\_docs&rev=1473428510](https://sermn.uab.cat/wiki/doku.php?id=informatica:hp_z240mt_biospec_docs&rev=1473428510)

Last update: **2016/09/09 15:41**

